




Hypertronics Male 9 Pin

Y- Probe



Pin Housing Rear View

(3 & 4) R1 = 1.3Ω k



REV	Date	N*	Drawn	Approv	Material :
					Part No. 0018602


Title P860YA Simed Baxter			
Scale Not To Scale		Dim in	mm
Date 04/11/03		Tol	± 0.2
Dwg No. SPY-860		Drawn	K.J. Rush





VIAMED Ltd.
15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT

Hypertronics Male 9 Pin


Y- Probe



Enlarged front view



Title P860YA
Simed
Baxter



Dim in mm

Tol ± 0.2

Drawn J.Nirwan

Scale Not To Scale

Date 10/05/02


Dwg No. SPY-860

VIAMED Ltd.
15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT


REV	Date	N°	Drawn	Approv	Material :
					Part No. 0018602

Hypertronics Male 9 Pin


Y- Probe



Enlarged front view




Hypertronics male 9 pin



R1 - 1.33 K OHM 1%, 1/8W RESISTOR

* INNER SHIELD AND OUTER SHIELD TO BE TWISTED TOGETHER.

Title P860YA
Simed
Baxter



Dim in mm

Tol ± 0.2

Drawn J.Nirwan

VIAMED
15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT



Ltd.

REV	Date	N°	Drawn	Approv	Material :
					Part No. 0018602

Scale	Not To Scale	Date	10/05/02	Dwg No.	SPY-860
-------	--------------	------	----------	---------	---------


Hypertronics Male 9 Pin

Y- Probe




Hypertronics Male 9 Pin

Y- Probe



Hypertronics male 9 pin



R1 = 1.33 k OHM 1%, 1/8W RESISTOR
 * R1 SHIELD AND OUTPUT SHIELD TO BE TWISTED TOGETHER