

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Valid from 01/01/2023

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8

Ceratherm 600-3 Radiant warmer with mounting systems



Pos.	RefNumber / GTIN	Description	Package	EURO/Unit

1	521A-60020-3-EX	Ceratherm 600-3 Radiant warmer		1'198.00
	7649994198015			



- Compact housing, color white RAL 9010
- · Ceramic heating element with 600 W
- · 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.

Product Information:

- Built-in two-stage LED lighting
- · Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- Without mounting bracket
- Weight: 4.3 kg

Further information on page 7 or in the

Data sheet Ceratherm 600-3 No. 15.249.xx.



Pos.

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland

Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8 Product Information:



EURO/Unit

Ceratherm 600-3 Radiant warmer with mounting systems

Description

Valid from 01/01/2023


Package

	GTIN			
2	521A-60023-3-EX	Ceratherm 600-3 Radiant warmer with fix wall mounting bracket (not pivotable)	1	1'335.00
	7649994198046			
2a	521-071204-EX	Fix wall mounting bracket (not pivotable), without radiant warmer	1	198.00
	76/000/108030			



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With fix wall mounting bracket (not pivotable).
 Dimensions as per drawing
- · Weight: 6.8 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





Pos.

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8



EURO/Unit

Description

Valid from 01/01/2023

	GIIN			
3	521A-60022-3-EX	Ceratherm 600-3 Radiant warmer with double-joint arm and wall mounting bracket	1	1'848.00
	7649994198060			
3a	521-093008300-EX	Double-joint arm and wall mounting bracket, without radiant warmer	1	875.00
	7649994198053			
3b	523-093.008.094-EX	Mounting plate for wall bracket (300x300x10mm, RAL 9010), optional accessories for lightweight construction	1	96.00
	7649994198138			



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.


Product Information:

Package

- Built-in two-stage LED lighting
- · Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm and wall mounting bracket. Dimensions as per drawing. Special lengths on request (*).
- Weight: 11.2 kg

Further information on page 7 or in the

Data sheet Ceratherm 600-3 No. 15.249.xx.





GTIN

Pos.

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8



EURO/Unit

Ceratherm 6	600-3 Radiant	warmer with	mounting	systems
-------------	---------------	-------------	----------	---------

Description

Valid from 01/01/2023

	<u> </u>				
4 (*)	521A-60025-3-EX	Ceratherm 600-3 Radiant warmer with double-joint arm and wall mounting bracket, with height adjustment (spring balanced arm)	1	2'179.00	
***************************************	7649994198121				
4a	523A-25-3-EX	Double-joint arm and wall mounting bracket, with height adjustment (spring balanced arm), without radiant warmer	1	1'445.00	
	7649994198077				
4b	523-093.008.094-EX	Mounting plate for wall bracket (300x300x10mm, RAL 9010), optional accessories for lightweight construction	1	96.00	
	7649994198138				


- Compact housing, color white RAL 9010
- · Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.

Product Information:

Package

- Built-in two-stage LED lighting
- · Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm, wall mounting bracket, and height adjustment. Dimensions as per drawing. Special lengths on request (*).
- Weight: 13.7 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





GTIN

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Valid from 01/01/2023

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Ceratherm 600-3 Radiant warmer with mounting systems

Discount scales page 8 Product Information:




Pos	RefNumber /	Description	Packago	EURO/Unit
Pos.	OTIN	Description	Package	EURU/UIII

5 (*)	521A-60024-3-EX	Ceratherm 600-3 Radiant warmer with double-joint arm and ceiling anchorage, with height adjustment (spring balanced	1	2'554.00
	7649994198091			
5a	523A-24-3-EX	Double-joint arm and ceiling anchorage, with height adjustment (spring balanced arm), without radiant warmer	1	2'016.00
	7649994198084			



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- · Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm, ceiling anchorage and height adjustment. Dimensions as per drawing.
 Special lengths on request (*).
- Weight: 30.0 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





7649994198107

GTIN

Pos.

Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8 Product Information:



EURO/Unit

Ceratherm 600-3 Radiant warmer with mounting systems

Description

Valid from 01/01/2023


Package

	• · · · ·			
6	521A-60035-3-EX	Ceratherm 600-3 Radiant warmer on mobile floor stand, with height adjustment	1	2'224.00
	7649994198114			
6a	523-90057150-EX	Mobile floor stand, with height adjustment, without radiant warmer	1	1'408.00



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- · Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With triple-wheel floor stand. Height adjustable. Two wheels with brake. Dimensions as per drawing.
- Weight: 37.5 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8 Product Information:



Ceratherm 600-3 Radiant warmer with mounting systems

Valid from 01/01/2023

Specifications Pos. 1 - 6 (**)

Model	CERATHERM
Туре	600-3
Article number	521A-60020-3
Conformity	CE mark (Directive 93/42/EEC)
Notified Body	TÜV SÜD Product Service GmbH, DE-80339 München (CE 0123)
IP protection class	IP20
Contamination level class	2
Electrical protection rating	I (with protective earth connection)
Medicinal product class	IIb (Directive 93/42/EEC, Annex IX, Rule 9)
Applied industrial standards:	- EN 60601-1:2006+A1:2013 - EN 60601-2-21:2009+A11:2011+A1:2016 - EN 60601-1-2:2015
In accordance with the Directive	2015/863/EU (RoHS-3)
Mains connection voltage	230 V +/-10%, 50-60 Hz
Power consumption, standby	0.03 VA
Power consumption, max.	630 VA
Fuses	2 × 3.15 AT microfuses, 5 x 20mm (in the unit socket) 1 x 800mAT microfuse, 5 x 20mm (internal on the powerboard)
External conditions during transport	- Temperature: −10 to +50 °C - Relative air humidity: 20 to 90% (no condensation)
External conditions during storage	- Temperature: -10 to +50 °C - Relative air humidity: 20 to 90% (no condensation)
External conditions during operation	- Temperature: +5 to +30 °C - Relative air humidity: 20 to 90% (no condensation)
Heating element	Ceramic, 600 W, wavelength 2-10µm (IR-B/C)
Lighting	LED, 1 x 3W / 700mA / 4.5V, colour temperature 4000°K, without built-in power source



Nufer Medical - let's find the best solutions

Nufer Medical AG

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

PRICE LIST for Distributors (not valid in Switzerland)

Discount scales page 8 Product Information:



Ceratherm 600-3 Radiant warmer with mounting systems

Valid from 01/01/2023

Alarms	 Acoustic and visual after 15 minutes, with automatic heat output reduction to 20% of max. output. Mains failure Technical malfunction (with error code)
Operating position	Horizontal, heat radiating vertically downwards
Assembly options	As per with accessories, chapter 16 / User Manual
Transport In original packaging	
Safety checks	Every 12 months as per service manual specifications, chapter 20
Disposal	In accordance with Directive 2002/96/EC. Get in touch with your allocated sales partner.
Dimensions	270 × 105 × 540 mm (W × H × L)
Weight	4.3 kg (without accessories)
Material, housing	ABS (temperature-resistant plastic)
Material, chassis	CNS (stainless steel)
Material, reflector	CNS (stainless steel)

Discount scales:

Units per order (Pos. 1 - 6)	1 – 4	List price
Units per order (Pos. 1 - 6)	5 – 9	List price - 3 % discount
Units per order (Pos. 1 - 6)	10 – 19	List price - 5 % discount
Units per order (Pos. 1 - 6)	20 – 49	List price - 8 % discount
Units per order (Pos. 1 - 6)	50 and more	List price - 10 % discount

^(*) Additional charge for special lengths at positions 4 (4a) and 5 (5a): € 199.00 per system.

^(**) NUFER MEDICAL reserves the right to make changes to the unit that may affect these specifications, at any time and without notice.