

# Training Feedback Form

Training Course Completed: <i>Temperature Probe Training</i>			
Date: <i>23/1/23</i>	Time/Length: <i>30/45mins</i>	Trainer: <i>Cabrin.</i>	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If so, what would you like this training to cover?			
Further comments:			
Name: <i>Aqib Majeed</i>			
Date: <i>23/1/23</i>			

\*Please delete as applicable



## Temperature Probes Basic Competency Questions

Please use your PCs as reference

1. What are the names of the connectors of the Temperature Probes we provide:

- Right Angled Jack - VSI 400 ✓
- Jack - GE Healthcare ✓
- Central - Philips/Hewlett Packard. ✓
- Side - Dräger/Siemens ✓

2. What is a brief description of a Temperature Probe and which department is it used in?

For Measurement of body Temperature and Core temperature in all ages and are usually used in Hospitals ✓

3. What area and which patient type can part number 0212016 be used on?

Throat / upper part behind nose, Paediatric Patient. oesophageal/rectal ✓

4. Can the Temperature Probes to be reused?

~~No~~ Yes ✓

5. Where can you find marketing materials?

can be found on stock page under 'i' icon. ✓

6. If you searched the online Temperature selection guide for the following, which part number would you get?

Application: Skin, Connector: HP/Philips, Patient: Neonatal

0212030 ✓

7. Do the temperature probes contain Latex?

No ✓

8. Find a part number and Euro € Distributor price for Temperature Probe Skin Contact Adult for Hewlett Packard (HP).

Part number: 0212032 ✓ Price: €54.00 ✓

9. On Intrastats, where can you find the Cross Reference?

under the Magnifying glass ✓

10. What date and part number is on the Temperature Probe leaflet?

Part number: 0290011. Date: 09/22 ✓