

7 May 1998

Part number	Compatible manufacturer	Dot colour	Pad seat colour	Pad colour	Cable			Comments
					Colour & material	Length		
						Feet (ft)	Metres (m)	
P856RA	Nellcor	Blue	Blue	White	White PVC	3	0.90	
P857RA	Spacelabs	White	Blue	White	White PVC	10	3.00	
P858RA	Nellcor	Blue	Blue	White	White PVC	12	3.65	
P859RA	Pace Tech	White	Blue	White	White PVC	12	3.65	
P860RA	Baxter	White	Blue	White	White PVC	12	3.65	
P861RA	BCI	Orange	Orange	White	White PVC	3	0.90	
P862RA	Criticare	Purple	Purple	White	White PVC	10	3.00	
P863RA	Datascope	White	Green	Black	White PVC	3	0.90	
P864RA	Datascope	White	Green	Black	White PVC	12	3.65	
P865RA	Sensormedics	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P867RA	Ohmeda	Green	Green	Black	White PVC	12	3.65	
P869RA	Criticare	Purple	Purple	White	White PVC	3	0.90	
P870RA	Critikon	White	White	White	White PVC	3	0.90	Discontinued, replaced by P865RA.
P871RA	Nonin	White	Tan	White	White PVC	3	0.90	
P872RA	Datex	Tan	Tan	White	White PVC	12	3.65	
P873RA	Datex	Tan	Tan	White	White PVC	3	0.90	
P875RA	Novametrix	Yellow	Yellow	Black	White PVC	10	3.00	
P876RA	Novametrix	Yellow	Yellow	Black	White PVC	10	3.00	
P877RA	Critikon	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P878RA	Invivo	White	Blue	White	White PVC	10	3.00	
P886RA	Kontron	?	?	White	White PVC	12	3.65	Not released yet due to unavailability of connectors.
P887RA	Spacelabs	Yellow	Yellow	Black	White PVC	12	3.65	Unissued part number.

7 May 1998

Part number	Connector				Supplied by MCI	Derivation	Comments
	Type	Number of contacts	Gender	Comments			
P856RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P857RA	Hypertronics	12	M	Pins 5,6,7,9 & 10 missing.		P858	
P858RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P859RA	Switchcraft	5	F			P858	
P860RA	Hypertronics	9	M	Pins 2 & 5 missing?		P858	
P861RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P862RA	Lemo	5	M		✓	Original	
P863RA	D	9	M	Pin 8 missing.		P867	
P864RA	DIN	8	M			P867	
P865RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P867RA	Hypertronics	9	M	Pins 3 & 6 missing.	✓	Original	
P869RA	D	9	M			P862	
P870RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P871RA	D (Amp)	9	M			P873?	
P872RA	DPK1	10	F	Pins 2,3 & 0 missing.	✓	Original	No connector supplied
P873RA	D	9	M	Pins 2 & 8 missing.	✓	Original	
P875RA	Hypertronics	7	M	Pin 7 missing.		P876	Wire cross over
P876RA	Hypertronics	7	M	Pin 7 missing.	✓	Original	
P877RA	D	9	M	Pins 7, 8 & 9 missing		Kit	Sensormedics LED's & cable.
P878RA	Lemo	7	M			P858	Query cable length
P886RA					✓	Original	
P887RA	Hypertronics	12	M	Pins 5,6,7,8 & 12 missing.		P876	

7 May 1998

Part number	Connection		Extension cable part number
	Monitor	Extension cable	
P856RA	✓	✓	Nellcor 958E8 or P958E4. Spacelabs P926E10
P857RA	✓		
P858RA	✓	✓	
P859RA	✓		
P860RA	✓		
P861RA		✓	
P862RA	✓		
P863RA		✓	P963E10
P864RA	✓		
P865RA		✓	
P867RA	✓		
P869RA		✓	P968E8 or P969E8
P870RA		✓	
P871RA		✓	P971E3 or P971E10
P872RA	✓		
P873RA		✓	P973E10
P875RA	✓	✓	Spacelabs P925E2
P876RA	✓	✓	Spacelabs P925E2
P877RA		✓	
P878RA	✓		
P886RA	✓		
P887RA	✓		