



~~ACT~~ VIAMED INVIVO HYPER-TRONICS PROBE


N.B. see previous pages for original probe.

LED face up.

PN: PDI-E833

FINGER
← PROBE.* LED 930nm = IR
660nm = R+This LED is the
same as the
Simer.

LED



ON

← Y

PROBE.

Y R

Photodiode R+

PN: PDB-C165




W B



CONNECTOR R+ 7PIN HYPERTRONICS (BLACK)

R1 = 39.2KΩ (1%)



- 1/ BLACK
- 2/ WHITE
- 3/ RED
- 4/ ORANGE
- 5/ MAIN SHIELD
- 6/ R1
- 7/ R1

LED PICKUP.

MAKE SURE LED IS WIRED CORRECTLY
ORANGE & RED FOR FINGER.
YELLOW & RED FOR Y.