



Physio Control Generic Probe



Nellco - Photo

Nellco - LED




ATT Steve

Here is Physio Control Diagram

Thanks Dave

PhysioControl

LED



Face Up.


PCB

Original

our cable

- | | |
|-----------------|--------|
| ① Yellow | Yellow |
| ② Red | Red |
| ③ Green | Orange |
| ④ Shield (main) | Shield |
| ⑤ Blue | Blue |


Sensor



Face Up.

- | |
|---------|
| ① White |
| ② Black |

Connector




(nos are not marked on connector)

Original


our cable

- | | |
|----------|----------------|
| ① Green | Orange |
| ② - | - |
| ③ Shield | Shield (inner) |
| ④ Red | Red |
| ⑤ Yellow | Yellow |
| ⑥ Black | Black |
| ⑦ Shield | Shield (main) |
| ⑧ Blue | Blue |
| ⑨ White | White |


①



④




Phys.cmx






- 1 Not Used
- 2 Green
- 3 Red
- 4 Coloured Shield MAIN
- 5 Blue
- 6 White
- 7 Black
- 8 B/W Shield INNER
- 9 Yellow

NB - resistor is PCB mounted at clip end.

Sensor Face Up



Led Face Up



Drawn By D.Lamb

Signed