


FIRST LOOK AT DISPOSABLE NEUCOR CABLES.

The red and orange LID wires have swapped place when compared to the P306 cable. The wiring including the shields is correct to the Neucor original. ~~7.5K Ω when measured~~
 The resistor is 7.5K Ω when measured with a multimeter.

Athens ext. cable for axi'hip.

EV



✓ 1 - 7

✓ 2 - 5

✓ 3 - NC

✓ 4 - B3 light is in small black connector

✓ 5 - B 8

✓ 6 - NC


✓ 7 - 6?

✓ ~~WXT~~ 8 - 6 Y 6

✓ 9 - Sh / W

✓ 0 - Bk 1

FV



1, Bk - (sen)

2, R - Com

3, R - Com

4, Sh - Sh (sen)

5, 9 - RED

6, Y - RES

7, 0 - R

8, Bk - RES

R X Y Bk

9 0 Bk S