



COMDEX MONITOR MODEL NO: MD-600P


OPEN the battery compartment and remove screen cover



Female conn. Pin




EXT. CABLE FOR EDENTEL monitor. REF SCP-10



Hand strip side.
Metal surround connects to the main shield.

Connector is 3M 10114
Reorder RS 299-4604

RV



Red wire side

PTC

EXPECTED COMDER FINGER PROBE

660nm (red)

905nm (IR)

MEASURED

650nm (red)

720nm (IR)

SCI MCI


260

715

SCI

661


712




RV MALE (RED)

- 1, WHITE
- 2, SHIELDS
- 3, BLACK
- 4, RED
- 5, -
- 6, GREEN
- 7, SHIELDS

CABLE LENGTH = 1.8m




COMDER MCI WIRE




- 1, WHITE
- 2, SHIELDS
- 3, YELLOW
- 4, RED
- 5, -
- 6, BLACK
- 7, SHIELDS





Redel Male 7 Pin




Finger Probe



Connector Rear View



Title P 85 4RA
Medlab/Com dek



Dim in

mm

Tol

± 0.2

Drawn

J.Nirwan

Scale

Not To Scale

Date

03/12/01


Dwg No.

SPF-854

Material:

Part No. 0018540

REV Date N° Drawn Approv



VIAMED Ltd.

15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT