



Redel Male 7 Pin

Finger Probe



Enlarged front view

Redel Male 7 Pin



* INNER SHIELD AND OUTER SHIELD TO BE TWISTED TOGETHER.

Title P855RA
BCI



Dim in mm
Tol ± 0.2

Scale Not To Scale



Drawn J.Nirwan

Date 01/05/02

Dwg No. SPF-855.1

Material :
Part No. 0018550

VIAMED Ltd.
15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT



Title P855RA BCI			
		Dim in	mm
		Tol	± 0.2
Scale	Not To Scale	Drawn	J.Nirwan
Date	01/05/02	Dwg No.	SPF-855



VIAMED Ltd.
15 Station Rd
Cross Hills, Keighley
West Yorkshire
BD20 7DT

Date	24/03/97	Type	BCI (EQUIVALENT E344)
bci.cmx		WORKS ON NELLCOR MACHINE & TEST BOX	

REAR VIEW

4ft

DISCARD
STRAIN RELIEF

On small porcelin clip style Sensor:

White goes onto small contact
on black coloured sensor surface

Black goes to larger contact on
reverse (Gold coloured)

SW 1

SW 4

Drawn By:	D.Lamb	Signed	
-----------	--------	--------	--