Service Repair Sheet SRS68558

Contact Name

Ashley Davies

Company/ Hospital Name

Thames Medical Ltd

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

01903522911

01903522911

00011982

ashley@thamesmedical.com

31/Aug/2023

Robert Connor

Viamed

Warranty



Notes 30/Aug/2023 Catherine Spence 30/Aug/2023 Catherine Spence Returning Unit repaired on SRS68547.

Connected Officted Off St S000547.

Capnograph S/N 81030156

Email states that when they received the unit, they checked it before sending out and that the unit did not power up. They are wondering if a ribbon issue . can we have a quick turn around

31/Aug/2023 Robert Connor

Received 1 x VM-2500-M-Vet s/n 81030156, with fitted mounting plate

Ready For quote

Repair Complete Signed

Caeen 5.9.23

SRN	Equipment	Stock Ref	Serial Number	Warranty
SRN35977	Capnograph	4410530	81030156	Yes (r)

4440000 × 1 foc SIN, SRS, SRN Repaired under warranty OPS × 1 foc

Thames Medical Ltd

8 A'Becket Gardens

Worthing, West Sussex BN13 2BW

www.thamesmedical.com

email

info@thamesmedical.com 01903 522911

tel

fax

01903 522922

Experts in patient monitoring and anaesthesia

Thames Medical

PURCHASE ORDER

order#

TOTAL

PO-01974

£0.00

date

30/08/2023

vendor address Viamed Ltd

15 Station Road

Cross Hills

Keighley, West Yorkshire

BD20 7DT

contact

Kate

phone fax

01535 634542

01535 635582

carrier

payment terms

Net 30

UPS

ITEM	DESCRIPTION	VENDOR PRODUCT CODE	QUANTITY	UNIT PRICE	SPECI AL TAX RATE	SUB-TOT
TML-V-VM250 0M1-RL2WR	VM2500M - Warranty Repair - Service Level 2 SN:81030156	Warranty Repair of VM2500M Mainstream Capnograph	1	£0.00	Taxable	£0.00
remarks	SRS68547			SUB-TOTAL		£0.00
	Upon return will not	power on under any circums	tance.	VAT		£0.00

Service Repair Sheet 68558

Contact Name

Ashley Davies

Company/ Hospital Name

Thames Medical Ltd

Direct Phone General Phone 01903522911

Opera Account

01903522911 00011982

Email

ashley@thamesmedical.com

Date Received

31/Aug/2023

Booked in By

Robert Connor

Repair

Ref

S/N

Equipment Type

Under Warranty

SRN35977

4410530 81030156

Capnograph

Complete - Repaired

Time: 0 Hour(s)

The Capnograph VM 2500-M Vet has been tested. The unit was powered up numerous times from mains PSU, AA batteries and a rechargeable battery. No issue was found with any power source.

A problem with the SpO2 sensor working intermittently arose. The SpO2 sensor connector has been rewired to correct this issue.