



Single Patient Use Oxygen Hoods

Designed for use with neonatal and paediatric patients, for the short-term, stable delivery of oxygen



These Oxygen Hoods consist of flat sheets of transparent vinyl with built-in support frames and a drape neck opening.

- See-through environment provides full patient visibility.
- Drape neck opening provides easy access to the patient.
- Single patient use means reduced risk of cross-infection.
- Flat packed for easy storage and transportation.

Standard Hood

The standard Oxygen Hood has two 22mm inlet ports, a full vinyl bottom, a drape neck opening and is available in a variety of sizes.

Other variants include:

- Partial bottom - remains stable whilst reducing heat build-up under the hood
- No-bottom - provides easier access to the patient

NICU


The NICU single patient use hoods are bottomless with a 75mm x 125mm draped access port on the rear panel, plus a drape to the front. Offers maximum patient access for free flow (blow-by) procedures or adjustments.

Double Oxygen Tent Set


The Double Tent Set is a combination of the Standard Hood with a specially designed outer tent with four 22mm ports for vacuum and aerosol tubing. Excess aerosolized medications are vacuumed from the space between the tents. The aerosol port of the inner tent is matched with the outer tent.

Low Birth Weight Oxygen Hood


For full-body oxygen/aerosol treatments, Low Birth Weight Oxygen Hoods are available in two sizes. These hoods can also be used in high-humidity/mist procedures. A variety of access port options include top access, front drape and rear access port. Low Birth Weight Oxygen Hoods are ideal for neonates on warming tables.




Standard Oxygen Hood




No Bottom Panel Oxygen Hood




Partial Bottom Panel Oxygen Hood



NICU Oxygen Hood



Double Oxygen Tent Set




Low Birth Weight Oxygen Hood




VIAMED


Single Patient Use Oxygen Hoods




Standard



Partial bottom




No bottom




NICU

Dimensions (mm) L x W x H	Standard	Partial bottom	No bottom	NICU
178 x 178 x 127			2310025	
178 x 178 x 178	2310000		2310026	2310075
216 x 216 x 216	2310001		2310027	2310076
254 x 254 x 203	2310003		2310028	2310077
254 x 254 x 254	2310004		2310029	2310078
305 x 305 x 254	2310006		2310031	2310080
305 x 305 x 305	2310007	2310050	2310032	2310081
457 x 457 x 457	2310010	2310052	2310035	
610 x 610 x 610	2310013	2310055		



Low Birth Weight Oxygen Hood

Dimensions (mm) L x W x H	Front drape 22mm inlet: 1 on LHS & 1 on RHS	Front Drape 22mm inlet: 1 on LHS & 1 on RHS Access port: Rear panel	Front drape 22mm inlet: 1 on LHS & 1 on RHS Access ports: Rear panel & top panel
254 x 406 x 127	2310100	2310101	2310103
279 x 470 x 178	2310104	2310105	2310107



Note: Excludes tubing, connectors and filters

Double Oxygen Tent Set

	Inner Tent Dimensions (mm) L x W x H			
		254 x 254 x 254	305 x 305 x 305	457 x 457 x 457
Outer Tent Dimensions (mm)	305 x 305 x 305	2310130		
L x W x H	457 x 457 x 457	2310129	2310127 2310128*	
	610 x 610 x 610		2310125	2310126

CE 0086

Specifications Subject to Change.

* No bottom panel

Distributed by:



VIAMED

Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582
Email: info@viamed.co.uk Website: www.viamed.co.uk

Part Number: 0392051
Date: 09/12



BS EN ISO 9001:2008
ISO 13485:2003

