pulse oximetry



MD300-C

Finger Pulse Oximeter

Ultra-compact, lightweight pulse oximeter with LED display

Low cost

Easy to use; 1-button operation

Automatic power-off when finger is removed

Economical pulse oximeter for general applications



Actual Size

The MD300-C can be used in a variety of applications, from use in hospitals, paramedic ambulances, clinics and dental surgeries, through to spot-checking before, during and after exercise.



MD300-C Finger Pulse Oximeter

Measuring Range

 SpO2
 70 - 100%

 Pulse Rate
 30 - 250 bpm

Accuracy

 $\begin{array}{lll} {\rm SpO}_2(70\text{ - }100\%) & \pm 2 \text{ digits} \\ {\rm SpO}_2 \text{ Resolution} & 1\% \\ {\rm Pulse \ Rate; \ beats \ per \ minute \ (30\text{ - }99\text{ bpm})} & \pm 2 \text{ bpm} \\ {\rm Pulse \ Rate; \ beats \ per \ minute \ (100\text{ - }250\text{ bpm})} & \pm 2\% \\ {\rm Pulse \ Rate \ Resolution} & 1 \text{ bpm} \end{array}$

Display

Features LED display

Data displayed Oxygen saturation, pulse rate, pulse indicator bar graph

Indicators Signal strength, battery status

Environmental Conditions

Operating temperature 5°C to 40°C Storage temperature -20°C to +55°C

Regulatory Approvals

Standards CE Class Ila

Assembly

Dimensions 56 mm (L) x 32 mm (W) x 30 mm (H)
Weight Approx. 50 g (including batteries)
Power Supply 2 x AAA/LR03/MN2400 1.5V alkaline batteries
Battery Life Approximately 16 hours of continuous operation.
Automatically powers off when finger is removed for longer than 8 seconds.

Part Numbers

MD300-C 2810010

EAN-13 (GTIN-13) Barcode: 5051826003046 Standard accessories included:

Lanyard

2 x AAA Batteries

Operator Manual

Optional AccessoriesSoft Carry Case, with belt loop, colour - black0022190Soft Carry Case, with belt loop, colour - orange0022191Soft Carry Case, with belt loop, colour - green0022192Soft Carry Case, with belt loop, colour - blue0022193

Clinical trials carried out in accordance with BS EN ISO 9919:2005 standards. Materials tested in accordance with ISO 10993-, -5, -10.







Optional carrying cases with belt loop and D-ring

€0123

Specifications Subject to Change

Distributed by:



NOTABLE

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk



