

CERATHERM 600-3

Radiant Heater




Shipping Dimensions / Packaging


Contents

1.1	WALL MOUNT, FIXED, NON-PIVOTING (OPTIONAL)	3
1.2	WALL MOUNT WITH DOUBLE-JOINT ARM (OPTIONAL)	4
1.3	WALL MOUNT WITH DOUBLE-JOINT ARM, LONG, WITH HEIGHT ADJUSTMENT (OPTIONAL)	5
1.4	CEILING MOUNT WITH DOUBLE-JOINT ARM, LONG, WITH HEIGHT ADJUSTMENT (OPTIONAL)...	6
1.5	MOBILE STAND, WITH HEIGHT ADJUSTMENT (OPTIONAL)	7



1.1 Wall mount, fixed, non-pivoting (optional)

	
Article number: 521A-60023-3-(CH/EX)	Non-pivotable. Distance from wall adjustable as per drawing.
Shipping Dimensions / Packaging:	
CERATHERM 600-3 radiant heater	1 box, 580 x 300 x 140 mm , weight 5.5 kg (Ceratherm and box)
Wall mount	1 box, 400 x 200 x 200 mm , weight 2.6 kg (Wall mount and box)
Whole system	2 boxes, weight 8.1 kg (Ceratherm, wall mount, boxes)
	(Pallet weight: 24 kg)
Maximum number of systems per Pallet	10



1.2 Wall mount with double-joint arm (optional)

	
Article number: 521A-60022-3-(OND)	Pivotal, 360°. No height adjustment. Arm length as per drawing.
Shipping Dimensions / Packaging:	
CERATHERM 600-3 radiant heater	1 box, 580 x 300 x 140 mm , weight 5.5 kg (Ceratherm and box)
Wall mount with double-joint arm	Assembled, 810 x 440 x 100 mm , weight 6.9 kg (delivered on pallet)
Whole system	1 pallet, 1200 x 800 x 300mm , weight 36.4 kg (Ceratherm, arm, pallet)
	(Pallet weight: 24 kg)
Maximum number of systems per Pallet	5


1.3 Wall mount with double-joint arm, long, with height adjustment (optional)

	
Article number: 521A-60025-3-(CH/EX)	Pivotal, 360°. With height adjustment. Arm length as per drawing.
Shipping Dimensions / Packaging:	
CERATHERM 600-3 radiant heater	1 box, 580 x 300 x 140 mm , weight 5.5 kg (Ceratherm and box)
Wall mount with double-joint arm	Assembled, 900 x 350 x 300 mm , weight 9.4 kg (delivered on pallet)
Whole system	1 pallet, 1200 x 800 x 450mm , weight 38.9 kg (Ceratherm, arm, pallet)
	(Pallet weight: 24 kg)
Maximum number of systems per Pallet	3

1.4 Ceiling mount with double-joint arm, long, with height adjustment (optional)

	
Article number: 521A-60024-3(CH/EX)	Pivotal, 360°. With height adjustment. Arm length as per drawing. Ceiling tube length customer-specific
Shipping Dimensions / Packaging:	
CERATHERM 600-3 radiant heater	1 box, 580 x 300 x 140 mm , weight 5.5 kg (Ceratherm and box)
Ceiling mount with double-joint arm	Partially assembled, 900 x 500 x 300 mm , weight 25.7 kg (delivered on pallet)
Whole system	1 pallet, 1200 x 800 x 450mm , weight 55.2 kg (Ceratherm, arm, pallet)
	(Pallet weight: 24 kg)
Maximum number of systems per Pallet	2

1.5 Mobile stand, with height adjustment (optional)

	
Article number: 521A-60035-3-(CH/EX)	With height adjustment. Triple-bar trolley base with two brakes. Dimensions as per drawing.
Shipping Dimensions / Packaging:	
CERATHERM 600-3 radiant heater	1 box, 580 x 300 x 140 mm , weight 5.5 kg (Ceratherm and box)
Telescopic tube	1 box, 1510 x 115 x 110 mm , weight 5.9 kg (Tube and box)
Chassis with wheels and bracket	1 box, 810 x 620 x 240 mm , weight 29.7 kg (Chassis, bracket and box)
Whole system	1 pallet, 1200 x 800 x 1700 mm , weight 65.1 kg (3 boxes and pallet, box with tube vertically on the pallet)
	(Pallet weight: 24 kg)
Maximum number of systems per Pallet	5

Notes:**NUFER MEDICAL AG**

Morgenstrasse 148
CH-3018 Bern / Switzerland

Tel. +41 (0)31 958 66 66
Fax +41 (0)31 951 46 73

info@nufer-medical.ch
www.nufer-medical.ch