


# New Products from Viamed...

## VM-2400 Haemoglobin Monitor

SpO<sub>2</sub>, non-invasive Haemoglobin and Methaemoglobin multi-parameter monitor

- Simple and intuitive operation
- Fast measurement (approximately 40 seconds)
- Detection of unphysiological MetHb
- No bio-contaminated waste
- Pain-free
- Long battery life (up to 300 measurements)
- Storage of up to 12,000 measurements



## SMARTsat®

High performance SpO<sub>2</sub> technology for OEMs

The very latest and innovative signal processing technologies and algorithms enable precise measurements, even under very difficult physiological conditions. SMARTsat technology has extremely low power consumption requirements, and is a real alternative for OEM customers whom require high-end pulse oximetry solutions for their applications.

## VM-2160 SMARTsat High Performance Pulse Oximeter

Designed with the most difficult conditions in mind

- For spot checking or constant monitoring of SpO<sub>2</sub> and pulse rate
- Now incorporating SMARTsat SpO<sub>2</sub> technology
- Clear and bright, colour TFT display
- Robust, ergonomic design with durable silicone protective cover
- Includes a SMARTsat soft silicone sensor\*
- Stores and restores up to 480 hours of trend data
- Adjustable high and low alarms
- Antibacterial keypad



# VM-2020

## NIBP and SpO<sub>2</sub> Monitor

Ergonomic and lightweight ambulatory monitor that combines the very latest in blood pressure and pulse oximetry monitoring technologies

Utilizing patient-friendly blood pressure measurement during inflation of the cuff, halves the measurement time to just 17 seconds.

- Choice of blood pressure measurement: either during inflation or standard deflation measurement
- Choice of day and night time operating modes
- Up to 400 measurements or 30 hours of operation
- SMARTsat high performance SpO<sub>2</sub> technology provides accurate and reliable measurements
- Soft and flexible silicone pulse oximetry sensor optimizes patient comfort



## NIBP OEM Modules

NIBP Modules with SpO<sub>2</sub> measurement for patient monitors


Available as either NIBP only or combined with a SMARTsat OEM III SpO<sub>2</sub> module.

- Choice of blood pressure measurement: either during inflation or standard deflation measurement
- Measurement artifact suppression
- SMARTsat high performance SpO<sub>2</sub> technology provides accurate and reliable measurements



## Gas Sampling Lines S Type

High performance and economical sampling lines that are not adversely affected by anaesthetic gases



## Gas Sampling Lines H Type

High performance sampling lines that remove humidity and are not adversely affected by anaesthetic gases



**Viamed Limited** · 15 Station Road · Cross Hills  
Keighley · West Yorkshire · BD20 7DT · United Kingdom  
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)

Part Number: 9990993

