
New Products from Viamed...

VM-2400 Haemoglobin Monitor

SpO₂, non-invasive Haemoglobin and Methaemoglobin multi-parameter monitor

- Simple and intuitive operation
- Fast measurement (approximately 40 seconds)
- Detection of unphysiological MetHb
- No bio-contaminated waste
- Pain-free
- Long battery life (up to 300 measurements)
- Storage of up to 12,000 measurements





SMARTsat®

High performance SpO₂ technology for OEMs

The very latest and innovative signal processing technologies and algorithms enable precise measurements, even under very difficult physiological conditions. SMARTsat technology has extremely low power consumption requirements, and is a real alternative for OEM customers whom require high-end pulse oximetry solutions for their applications.

VM-2160 SMARTsat High Performance Pulse Oximeter

Designed with the most difficult conditions in mind

- For spot checking or constant monitoring of SpO₂ and pulse rate
- Now incorporating SMARTsat SpO₂ technology
- Clear and bright, colour TFT display
- Robust, ergonomic design with durable silicone protective cover
- Includes a SMARTsat soft silicone sensor*
- Stores and restores up to 480 hours of trend data
- Adjustable high and low alarms
- Antibacterial keypad







VM-2020

NIBP and SpO₂ Monitor

Ergonomic and lightweight ambulatory monitor that combines the very latest in blood pressure and pulse oximetry monitoring technologies

Utilizing patient-friendly blood pressure measurement during inflation of the cuff, halves the measurement time to just 17 seconds.

- Choice of blood pressure measurement: either during inflation or standard deflation measurement
- Choice of day and night time operating modes
- Up to 400 measurements or 30 hours of operation
- SMARTsat high performance SpO₂ technology provides accurate and reliable measurements
- Soft and flexible silicone pulse oximetry sensor optimizes patient comfort





NIBP OEM Modules


NIBP Modules with SpO₂ measurement for patient monitors

Available as either NIBP only or combined with a SMARTsat OEM III SpO₂ module.

- Choice of blood pressure measurement: either during inflation or standard deflation measurement
- Measurement artifact suppression
- SMARTsat high performance SpO₂ technology provides accurate and reliable measurements

Gas Sampling Lines S Type

High performance and economical sampling lines that are not adversely affected by anaesthetic gases



Gas Sampling Lines H Type

High performance sampling lines that remove humidity and are not adversely affected by anaesthetic gases



Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk



Part Number: 9990993

