


capnography mounting options

for use with  **VIAMED** VM-2500 series of Capnographs



- 4410500** VM-2500-M - Central European
- 4410501** VM-2500-M - Asia
- 4410530** VM-2500-M - Veterinary




- 4410520** VM-2500-S - Central European
- 4410521** VM-2500-S - Asia
- 4410540** VM-2500-S - Veterinary



- 4410550** VM-2500-MG - Central European
- 4410560** VM-2500-MG - Veterinary



Universal Mounting Plate
for VM-2500 Range
4420600 Rev. 1
4420601 Rev. 2



0121180 V Wall Mount



0121181 V Mount Pole Clamp - Small
Clamping Width: 14 -25mm



0121182 V Mount Pole Clamp - Large
Clamping Width: 16 - 40mm



0121184 V Mount Rail Clamp
*Attaches to most types of
30mm (H) x 10mm (D) Medirail*



0121200 Universal Mounting Clamp
*Attaches to round or flat objects up
to 3.8 cm thick. Facilitates positioning
of the device at any angle*



Part Number: 0092201
Date: 05/13