Posey Products (Manufactured by TIDI Products LLC, USA)

Posey Pulse Oximetry Sensor Wrap

Posey sensor wraps are manufactured from soft foam and are used to hold pulse oximetry sensors in place.

They feature a hook and loop fastener that removes the need for adhesive tape, which can cause damage to fragile skin.




They are used with reusable wrap sensors and Y-sensors.



Posey wraps are not required with disposable sensors, as those feature self-adhesive tape as part of the sensor.

Posey sensor wraps can be used on various monitoring sites, depending upon patient size:

- Index finger (adult and paediatric patients)
- Foot (infant and small paediatric patients)
- Big toe (neonates)



The manufacturer (TIDI Products) recommends changing the pulse oximetry monitoring site every 2 hours, or in line with the recommendations of the pulse oximetry manufacturer or the local institution's protocols - this is to prevent skin damage and potential burns.

- Single patient use
- Supplied in boxes of 12
- Available in cases of 48 or 96 boxes of 12

Posey I.D. Bracelets

Posey I.D. Bracelets are soft, foam wraps that are used to provide identification of the patient without to use of adhesive tape.

Temporarily replaces permanent I.D. bracelet when permanent I.D. interferes with I.V. site.



Used in neonatal applications, available in 2 sizes:

- **Newborn** (14 cm x 3 cm)
- **Premie** (7 cm x 3 cm)

Posey I.D. Bracelets are:

- Easily moved and can be positioned on the wrist or ankle using hook and loop fasteners
- Single patient use
- Supplied in boxes of 12
- Newborn size is available in cases of 48 or 96 boxes of 12

Posey SecureSleeve

Multi-purpose external arm splint designed to help reduce patients' arm movements without the use of rigid splints



Used to hold the arm in place whilst an I.V. line is in place, or for the prevention and release of limb **contractures***.

*Contractures are a permanent tightening of the muscles and tendons that causes joints to shorten and stiffen.

- Supplied singly
- Available in the following sizes: Infant, Infant (large), Paediatric (2-5 years), Paediatric (2-7 years)
- · Can be used on hand, elbow, wrist and knee
- X-Ray translucent; can be left in place during X-ray procedures
- Reusable: see Instructions For Use for cleaning guidelines
- Optional SecureSleeve Clips are available to attach to the patient's sleeve and help prevent movement of the sleeve

Posey Elbow Splint

For use with patients that require restricted arm movement at the elbows to prevent self-injury (e.g. impairment of wound healing)



- Supplied as a pair / pack of 2
- Wraps around the arm at the elbow
- Uses hook and loop fastener to close
- Available in the following sizes: Infant (up to 2 years), Paediatric (2-7 years)
- Reusable: see Instructions For Use for cleaning guidelines
- · Contains stiffeners that can be removed for laundering

Posey I.V. Armboards

Flexible padded board, covered with soft vinyl for patient comfort.

Helps prevent arms from bending during I.V. procedures.



- Supplied in packs of 12
- Available in infant & paediatric sizes
- · Single patient use, but can be cleaned for re-use on the same patient
- Held in place by tape

Posey Foam Trach Ties

Recommended for use with patients that have **tracheostomy*** tubes or require oxygen cannula support.

*A tracheostomy is an opening created at the front of the neck so that a tube can be inserted into the trachea to assist breathing.







- Designed to help prevent extubation during severe coughing spasms
- Available in the following sizes: small, medium, large and extra large
- One-piece collar helps secure most tracheostomy tubes
- · Adjustable sizing for a comfortable fit
- Single patient use
- Supplied in packs of 12

Posey Secure Trach Ties

Recommended for use with patients that have tracheostomy tubes.

Features an extra tie to ensures that the tracheostomy tube is held in place more securely.



- Adult size: to fit necks 33 cm 48 cm circumference
- · Adjustable sizing for a comfortable fit
- Single patient use
- Supplied in packs of 12

Latex

All of the products that Viamed supplies are latex-free.

Where to find additional information

- Viamed website
- Product leaflets linked to stock pages
- FAQs on the stock pages
- Memos on the stock pages
- Instructions for use

Important Note: Some Posey products have been discontinued and are only available whilst stocks last.

Please check the stock pages in Intrastats for the latest information.