



VIAMED



Tom Thumb Test Sheet

Equipment

Description	Tom Thumb Vent.	Model	TT480
Date	23-02-00	S/N	9001 VM9001

Test Results

Test Equipment	Pressure Gauge Test	Pressure 10CmH20	Pressure 50 CmH20	P/F
P200L	Readings +/- 2			P
	Adjustable Valve @ 5 Lpm	Minimum 110 (<8 Cm H20) 9.8	Maximum SLTR 3CTR 46 (>43 & <47 Cm H20) 45	P
	Safety Valve @ 5 Lpm		Maximum (>42 & < 60 Cm H20) 50	P
Visual Checks				
Check	Gauge cannot be removed without tools.			P
Check	Gauge appears Straight.			P
Adjustable Valve set to Zero				P

Signed : _____