

VIAMED

Suppliers of Electromedical Equipment

2017



Ingleborough Falls, Ingletton

January

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	25
F		6	13	20	26
S		7	14	21	28
S	1	8	15	22	29

February

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	
T	2	9	16	23	
F	3	10	17	24	
S	4	11	18	25	
S	5	12	19	26	

March

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	31
S	4	11	18	25	
S	5	12	19	26	

April

M		3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

May

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

June

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	
S	4	11	18	25	

July

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

August

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

September

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	

October

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	25
F		6	13	20	26
S		7	14	21	28
S	1	8	15	22	29

November

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	
S	4	11	18	25	
S	5	12	19	26	

December

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31



VIAMED

15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk



VM-2500 Series Capnographs

To address the growing emphasis on the use of Capnography during resuscitation and intubation,

Viamed offer the VM-2500 range of high performance Capnographs

Available as Mainstream and Sidestream units; along with a range of VersaStream one-piece sampling lines (sidestream only)



Full scale accuracy within 10 seconds

VersaStream® CO₂ Sampling Lines

The alternative to Nafion® & Water Traps


Innovative sidestream monitoring solutions for use with a range of manufacturers' capnographs, including Philips, GE, Respirationics, Dräger and Viamed...

A one-piece sampling solution with integrated bacterial filter with active & complete removal of moisture and water, providing increased performance and reducing cost.


To find out more visit our Cross Reference
www.viamed.co.uk/linesxref



Nasal Sampling Line



Nasal Cannula



Nasal/Oral Cannula



Airway Adapter Sampling Line



Respironics/Philips



Luer Lock Male



Viamed VM-2500 series



VIAMED

Tel: 01535 634 542

Specifications Subject to Change
Part Number: 0196000

