

Compatibility Statement

To who is concerned

We, Medlinket declared that Disposable Non-invasive EEG sensor are compatible with all these monitor model# which are BIS module listed here (adapter cable required):




Model/Type reference: B-BIS-4A, B-BIS-4P






Issued BY: Shenzhen Med-link Electronics Tech Co., Ltd









Date: April 22, 2024.











End of this report






Compatible List of monitor brands and models




NO.	Compatible Brand	Model# (picture)
1	GE Healthcare	<div><p>CARESCAPE series: B450,B650,B850 etc.</p></div> <div><p>DASH series: B20,B40,B105, B125,B155 etc monitor.</p></div> <div><p>Datex Ohmeda Monitors S/5.</p></div>

NO.	Compatible Brand	Model# (picture)
2	Medtronic-Covidien	 <p data-bbox="592 645 1297 680">Covidien BIS VISTA/BIS VIEW Monitoring System.</p>
3	Mindray	 <p data-bbox="592 1059 983 1095">BeneVision N series monitor.</p> <p data-bbox="1070 1059 1442 1095">BeneView T series monitor.</p>
4	Comen	 <p data-bbox="592 1435 708 1471">K series.</p>  <p data-bbox="592 1753 868 1789">C series etc monitor.</p>  <p data-bbox="592 2029 735 2065">NC series.</p>

NO.	Compatible Brand	Model# (picture)	
5	Edan	 <p data-bbox="600 539 919 577">IX series (IX15/12/10).</p>	 <p data-bbox="1038 539 1469 577">Elite V series (V8/5/5) monitor.</p>
6	Drager	 <p data-bbox="600 943 879 981">Infinity Kappa series.</p>	 <p data-bbox="1038 943 1342 981">Infinity Gamma series.</p>
		 <p data-bbox="600 1312 935 1350">Delta series,Vista series.</p>	 <p data-bbox="1038 1312 1374 1350">Vista 120 series monitor.</p>
7	Nihon Kohden	 <p data-bbox="600 1895 919 1973">BSM-6000C series (6301C/6501C/6701C).</p>	 <p data-bbox="999 1895 1477 1973">BSM-1700 series (BSM-1733/BSM-1753/BSM-1763).</p>

NO.	Compatible Brand	Model# (picture)	
8	Philips		
		<p data-bbox="598 521 874 589">Intellivue MP Series Patient Monitor</p>  <p data-bbox="598 851 898 918">Philips M1204A V24C patient monitor</p>	<p data-bbox="1038 521 1313 589">Intellivue MX Series Patient Monitor</p>  <p data-bbox="1038 862 1452 896">M1205A V26C patient monitor</p>
9	Biolight	 <p data-bbox="598 1153 845 1187">A series(A8/6/5/3)</p>	 <p data-bbox="1038 1153 1268 1187">Q series(Q7/5/3)</p>
		 <p data-bbox="598 1429 917 1462">P series(P22/18/15/12)</p>	
10	Spacelabs	 <p data-bbox="598 1769 869 1803">Ultraview SL Series</p>	 <p data-bbox="1038 1769 1260 1803">90369 Ultraview</p>
		 <p data-bbox="598 2049 826 2083">91393 Xprezzon</p>	

NO.	Compatible Brand	Model# (picture)
11	BPL Medical	 <p data-bbox="596 510 1161 548">Modular multipara monitor ExcelSign E12</p>
12	Advanced	 <p data-bbox="596 913 890 952">Advanced PM-2000M</p>
13	Mennen medical	 <p data-bbox="596 1312 1098 1350">Menntor X7 Modular Patient Monitor</p>
14	Olidef Medical	 <p data-bbox="596 1680 1161 1718">Multi-parameter ECG monitor LOGICARE</p>
15	SINOHERO touching for life	 <p data-bbox="596 2045 962 2083">AcuitSign M series(M8/6/3)</p>

NO.	Compatible Brand	Model# (picture)
16	Innomed medical	
		P series(Anyview P18/15) Innocare T12 Plus/T17 Plus
17		
		Compact multi-parameter monitor SIGNSCOPE 15/12

*Disclaimer: All the registered trademarks, product names, models, etc. shown in the catalogue are owned by the original holder or the original manufacturer. This article is only used to explain the compatibility of the MED-LINKET products, and nothing else! All the above information is for reference only, and should not be used as a working guide for medical institutions or related units. Otherwise, any consequences will be irrelevant to the company.

Shenzhen Med-link Electronics Tech Co., Ltd

Address: 2th, 4th and 5th Floor, Building Two, Hualian Industrial Zone,
Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel: 86-0755-61120085 Fax: 86-0755-61120055
E-mail: marketing@med-linket.com
Website: www.med-linket.com

