

SpiroTrue H - OEM disposable single patient adult flow sensor

Technical data:

- Sensor measuring principle: differential pressure measurement
- Sensor reproducibility: $\leq 10\%$ at flow range 0 to 100 L/min
- Sensor designed acc 93/42/EEC (Anex I and Anex II)
- Sensor housing material: ABS
- Sensor weight: $\leq 9,5\text{g}$
- Sensor dead space: $\leq 10\text{mL}$
- Sensor connectors: 15M/15M - ISO 5356
- Sensor conforms with EN ISO 10993-5
- Sensor is none sterile packed

Storage conditions:

- Recommended storage temperature : -20°C to 50°C
- Humidity (r.H.) : 5% to 95%
- Recommended period of storage: ≤ 3 years
(by none opened packaging and keep away from sunlight)

Operating conditions:

- Recommended operating temperature : 10°C to 40°C
- Relative humidity $\leq 95\%$, none condensing


Notes:

- Before using of the sensor, consult the appropriate ventilator user instructions for more detailed information
- SpiroTrue H is compatible with Hamilton single use flow sensor PN: 279331

Product variants:

P/N: 1030132000 - Box of 6pcs: SpiroTrue H with 2m double tube (incl. hook and adaptor)

P/N: 4030131003 - SpiroTrue H with 2m double tube

										
					Date	Name	SpiroTrue H			
				Drawn by	12.07.2013	N.Radev				
				Checked by	15.07.2013	B.Lindner				
				Approved by	15.07.2013	B.Lindner				
							Data sheet -22102015		1	
									A3	
Status	Changes	Date	Name							