


patient monitoring

Pippa Breathing Monitor

- Instantaneous or averaged breathing rate.
- Easy to use and accurate.
- Pre-set alarm limits.
- Dedicated masks with integrated sensor; universal sensor available for use with other masks.
- Alternative to human observation or expensive multi-parameter patient monitors.



Utilising unique patented PEP sensors the Pippa breathing monitor reliably and accurately monitors patients' breathing rate.

 VIAMED

Pippa Breathing Monitor

Applications for use:

- General patient monitoring
- Oxygen therapy
- Epidurals
- Continuous spinals
- PCA
- Nebulised medication

Areas for use include:

- Accident and Emergency
- Paramedic Ambulances
- Acute Medical Wards
- Respiratory Units
- Obstetric Units
- High Dependency Units
- Postoperative Wards
- Recovery Wards
- Operating Theatres

Display

Breathing rate: 2 digits

Resolution: 1 breath per minute

Accuracy

+/- 1 breath per minute

Modes of operation

Instantaneous breathing rate

Averaged breathing rate

Display indications

AVERAGE BPM

INSTANTANEOUS BPM

BATTERY LOW

Alarm activation

Audio alarm triggered when two consecutive breaths are monitored

- i) below 8 b.p.m., 20 seconds between breaths
- ii) above 30 b.p.m., 4.8 seconds between breaths
- or
- iii) when battery is low

Alarm deactivation

Alarm is deactivated when:

- i) rate returns within normal operational parameters, between 8-30 b.p.m.
- ii) Rate Select button is pressed during an alarm

Dimensions

14.0 cm (l) x 6.4 cm (w) x 3.2 cm (h)

Weight

130g

Batteries

2 x AA, 1.5 V



0086

Patent Pending


The Pippa breathing monitor is an hand held device, but can be easily attached to a pocket or fastened onto a belt by using the integral securing clip.

Pippa monitoring kit includes:

Qty	Part number	Description
1	1610100	Pippa breathing monitor
1	1610200	Patient lead, monitor-sensor
1	1620010	Carrying case
1pk	1655000	Battery pack, AA, 1.5V. Pack of 2
1pk	1650010	Sensor, (universal). Pack of 5 x 1650009
1pk	1650050	Mask (Flexicare) with integrated sensor. Pack of 5 x 1650049
1	1690000	Instruction manual



Mask with integrated PEP sensor.



Universal and integrated PEP sensors.

Specifications subject to change.

Please contact Viamed for further specifications and information.



Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634542/636757 Fax: +44 (0)1535 635582
Email: info@viamed.co.uk Website: www.viamed.co.uk



BS EN ISO 9001
BS EN 46001
BS EN ISO 13485