

**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

### PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020

### Ceratherm 600-3 Radiant warmer with mounting systems

| Pos. | Item no         | Description                    | Package | EURO/Unit |
|------|-----------------|--------------------------------|---------|-----------|
|      | T               |                                |         |           |
| 1    | 521A-60020-3-EX | Ceratherm 600-3 Radiant warmer | 1       | 1'071.00  |



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- Without mounting bracket
- · Weight: 4.3 kg

Further information on page 7 or in the **Data sheet Ceratherm 600-3 No. 15.249.xx**.



**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

### PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020


#### Ceratherm 600-3 Radiant warmer with mounting systems

| Pos. | Item no         | Description   | Package | EURO/Unit |
|------|-----------------|---|---------|-----------|
|      |                 | Ceratherm 600-3 Radiant warmer with fix wall mounting             |         |           |
| 2    | 521A-60023-3-EX | bracket (not pivotable)   | 1       | 1'194.00  |
| 2a   | 521-071204-EX   | Fix wall mounting bracket (not pivotable), without radiant warmer | 1       | 170.00    |



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- · Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With fix wall mounting bracket (not pivotable).
   Dimensions as per drawing
- · Weight: 6.8 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

## PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020

### Ceratherm 600-3 Radiant warmer with mounting systems

| - | Pos.  | Item no          | Description  | Package | EURO/Unit |
|---|-------|------------------|--|---------|-----------|
|   | 3 (*) | 521A-60022-3-EX  | Ceratherm 600-3 Radiant warmer with double-joint arm and wall mounting bracket | 1       | 1'653.00  |
|   | 3a    | 521-093008300-EX | Double-joint arm and wall mounting bracket, without radiant warmer             | 1       | 754.00    |



- Compact housing, color white RAL 9010
- · Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm and wall mounting bracket. Dimensions as per drawing. Special lengths on request (\*).
- Weight: 11.2 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.

350 350 350 State of all left of the left



**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

### PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020


### Ceratherm 600-3 Radiant warmer with mounting systems

| Pos. Item no Description |                 | Package  | EURO/Unit |          |
|--------------------------|-----------------|--|-----------|----------|
| 4 (*)                    | 521A-60025-3-EX | Ceratherm 600-3 Radiant warmer with double-joint arm and wall mounting bracket, with height adjustment (spring balanced arm) | 1         | 1'949.00 |
| 4a                       | 523A-25-3-EX    | Double-joint arm and wall mounting bracket, with height adjustment (spring balanced arm), without radiant warmer             | 1         | 1'244.00 |



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- · 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm, wall mounting bracket, and height adjustment. Dimensions as per drawing. Special lengths on request (\*).
- Weight: 13.7 kg

Further information on page 7 or in the **Data sheet Ceratherm 600-3 No. 15.249.xx.** 





**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

## PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020

### Ceratherm 600-3 Radiant warmer with mounting systems


| Pos.  | Item no         | Description  | Package | EURO/Unit |
|-------|-----------------|--|---------|-----------|
| 5 (*) | 521A-60024-3-EX | Ceratherm 600-3 Radiant warmer with double-joint arm and ceiling anchorage, with height adjustment (spring balanced arm) | 1       | 2'285.00  |
| 5a    | 523A-24-3-EX    | Double-joint arm and ceiling anchorage, with height adjustment (spring balanced arm), without radiant warmer             | 1       | 1'736.00  |



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- · 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- · Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With double-joint arm, ceiling anchorage and height adjustment. Dimensions as per drawing.
   Special lengths on request (\*).
- Weight: 30.0 kg

Further information on page 7 or in the

Data sheet Ceratherm 600-3 No. 15.249.xx.



Price list **valid from 01.01.2020.** Replaces all previously issued price lists. All prices in EURO and do not include VAT. Handling and packaging costs included. Incoterm EXW 3018 Bern. Terms and conditions 10.007.K. Design and prices are subject to change without notice.

Terminal board for connecting SPR to electroinstallation



**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

### PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020


### Ceratherm 600-3 Radiant warmer with mounting systems

| Pos. | Pos. Item no Description |  | Package | EURO/Unit |
|------|--------------------------|--|---------|-----------|
|      |                          |  |         |           |
| 6    | 521A-60035-3-EX          | Ceratherm 600-3 Radiant warmer on mobile floor stand, with height adjustment | 1       | 1'989.00  |
| 6a   | 523-90057150-EX          | Mobile floor stand, with height adjustment, without radiant warmer           | 1       | 1'213.00  |



- Compact housing, color white RAL 9010
- · Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- With triple-wheel floor stand. Height adjustable. Two wheels with brake. Dimensions as per drawing.
- Weight: 37.5 kg

Further information on page 7 or in the Data sheet Ceratherm 600-3 No. 15.249.xx.





#### **Nufer Medical AG**

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

# PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Valid from 01/01/2020

## Ceratherm 600-3 Radiant warmer with mounting systems

Pos. Item no Description Package EURO/Unit

#### Specifications Pos. 1 - 6 (\*\*)

| Model                                | CERATHERM  |  |
|--------------------------------------|--|--|
| Туре                                 | 600-3  |  |
| Article number                       | 521A-60020-3   |  |
| Conformity                           | CE mark (Directive 93/42/EEC, Annex II)  |  |
| Notified Body                        | TÜV SÜD Product Service GmbH, DE-80339 München (CE 0123)   |  |
| IP protection class                  | IP20   |  |
| Contamination level class            | 2  |  |
| Electrical protection rating         | I (with protective earth connection)   |  |
| Medicinal product class              | IIb (Directive 93/42/EEC, Annex IX, Rule 9)  |  |
| Applied industrial standards:        | - EN 60601-1:2006 + A1:2013<br>- EN 60601-2-21:2009 + A11:2011 + A1:2016<br>- EN 60601-1-2:2015                      |  |
| In accordance with the Directive     | 2011/65/EU (RoHS II)   |  |
| Mains connection voltage             | 230 V +/-10%, 50-60 Hz   |  |
| Power consumption, standby           | 0.03 VA  |  |
| Power consumption, max.              | 630 VA   |  |
| Fuses                                | 2 × 3.15 AT microfuses, 5 x 20mm (in the unit socket)<br>1 x 800mAT microfuse, 5 x 20mm (internal on the powerboard) |  |
| External conditions during transport | - Temperature: −10 to +50 °C<br>- Relative air humidity: 20 to 90% (no condensation)                                 |  |
| External conditions during storage   | - Temperature: -10 to +50 °C<br>- Relative air humidity: 20 to 90% (no condensation)                                 |  |
| External conditions during operation | - Temperature: +5 to +30 °C<br>- Relative air humidity: 20 to 90% (no condensation)                                  |  |
| Heating element                      | Ceramic, 600 W, wavelength 2-10µm (IR-B/C)   |  |
| Lighting                             | LED, 1 x 3W / 700mA / 4.5V, colour temperature 4000°K, without built-in power source                                 |  |



**Nufer Medical AG** 

Morgenstrasse 148 PO Box 673 CH-3018 Bern Switzerland Phone E-Mail Internet VAT-N° +41 (0)31 958 66 66 info@nufer-medical.ch www.nufer-medical.ch CHE-106.962.174 MWST

Certified according to ISO 13485

## PRICE LIST FOR DEALERS (not valid in Switzerland)

Discount gradation see last page

Pos. Item no

Valid from 01/01/2020

**EURO/Unit** 

Package

#### Ceratherm 600-3 Radiant warmer with mounting systems

**Description** 

| Alarms              | <ul> <li>- Acoustic and visual after 15 minutes, with automatic heat output reduction to 20% of max. output.</li> <li>- Mains failure</li> <li>- Technical malfunction (with error code)</li> </ul> |
|---------------------|---|
| Operating position  | Horizontal, heat radiating vertically downwards   |
| Assembly options    | As per with accessories, chapter 16 / User Manual   |
| Transport           | In original packaging   |
| Safety checks       | Every 12 months as per service manual specifications, chapter 20  |
| Disposal            | In accordance with Directive 2002/96/EC. Get in touch with your allocated sales partner.  |
| Dimensions          | 270 × 105 × 540 mm (W × H × L)  |
| Weight              | 4.3 kg (without accessories)  |
| Material, housing   | ABS (temperature-resistant plastic)   |
| Material, chassis   | CNS (stainless steel)   |
| Material, reflector | CNS (stainless steel)   |

#### **Discount Graduation:**

| Units per order (Pos. 1 - 6) | 1 – 4       | List price                 |
|------------------------------|-------------|----------------------------|
| Units per order (Pos. 1 - 6) | 5 – 9       | List price - 3 % discount  |
| Units per order (Pos. 1 - 6) | 10 – 19     | List price - 5 % discount  |
| Units per order (Pos. 1 - 6) | 20 – 49     | List price - 8 % discount  |
| Units per order (Pos. 1 - 6) | 50 and more | List price - 10 % discount |

<sup>(\*)</sup> Additional charge for special lengths at positions 3 (3a), 4 (4a) and 5 (5a): € 163.00 per system.

<sup>(\*\*)</sup> NUFER MEDICAL reserves the right to make changes to the unit that may affect these specifications, at any time and without notice.