

12th November 2021

To Whom It May Concern

Ref. Clinell Universal Wipes
Clinell Sporidical Wipes

The TOF3D Neuromuscular Transmission Monitor is CE certified and details for cleaning the device can be found in the associated instructions for use.


Sales to the NHS require devices to be cleaned and disinfected with specific brands of products, which are readily available through NHS Supplies. Therefore, we have submitted several relevant types and brands of cleaning products to MIPM, the manufacturer of the TOF3D.

We have been advised that the Clinell products have successfully passed testing by MIPM and the associated report will be supplied to us in the near future. In due course the technical files and the TOF3D instructions for use will also be updated accordingly.

In the meantime should you have any queries or require further information, please do not hesitate to contact me.

Yours faithfully

Yours faithfully



Steve Nixon
Director – Viamed Ltd.