


RBSUS Flowmeter



Flowmeter — 0-5 Lbs/min
 No hose or hose coupling
 Blow off with loose profile
 Blow off safety
 No unnecessary holes

Vertical holes on Rear

Block Holes long

less than front

DRFOI



VIAMED



Drager Tom Thumb Resuscitation Unit Pt.No: DTT490

The basis of this configuration is a standard version of the Viamed TT490 Tom Thumb.

It comprises:-

The main block ,

0-5 ltrs/min Flowmeter no hose,



Safety valve

Adjustable valve

15 mm Male outlet.

Addition holes will be blanked off with blanking plates

Holes in the reverse for mounting Drager Bracket.



Viamed Limited, 15 Station Road, Cross Hills,
Keighley, West Yorkshire BD20 7DT
Tel: +44 (0) 1535 634542/636757 Fax: +44 (0) 1535 635582
Registration No. 1291765 in England

DTT490