Invoice Address Industrial Pipefreezing Services Ltd 8 The Coles Shop Merton Abbey Mills Watermill Way London **SW19 2RD**

Delivery Address Industrial Pipefreezing Services Ltd 8 The Coles Shop Merton Abbey Mills Watermill Way London SW19 2RD

Supplier Viamed Ltd 15 Station Road Cross Hills Cross Hills
Keighley, West Yorkshire
BD20 7DT, United Kingdom
Tel: +44 (0) 1535 634542
Fax: +44 (0) 1535 635582
Email: info@viamed.co.uk
VAT Reg No: GB287389593
Company Reg No: 01291765
Eori No: GB287389593000

Steve Turner Contact Name 02085439390 Contact Tel 00009085 Account Customer Reference 19042273ZT4 Date 20 Apr 2022

Tracking Number 1Z9W96386878217347

Invoice RVM136553-1 Paid

CPT Carriage Paid To United Kingdom * Incoterms® 2020

Delivery Reference DVM136553-1 Contact zoey.teal@viamed.co.uk

Item Reference	Description	Quantity	£ Unit	£ Unit Vat	£ Total
0131202 Tariff 90181990-00 CoO United Kingdom	Oxygen Sensor Cable Black - Viamed version Coiled cable with right angled connectors	1	20.00	4.00	24.00
0110017 Tariff 9019209000 CoO United States	Teledyne Oxygen Sensor R-17MED S/N:696887-696889	3	38.00	7.60	136.80
0100172 Tariff 9018199000	Mounting stand for AX300 & MX300 Excludes rubber feet (0131205 x 3 required per	1 unit)	20.00	4.00	24.00
0131205 Tariff 90181990-00 CoO U.S.A. PPUPS1	Suction cup foot - AX300 & MX300 base	2	5.00	1.00	12.00
	UPS Courier Delivery - Standard 1z9w96386878217347		0.00	0.00	0.00

Total Net: £ 164.00 Total Vat: £ 32.80 Total: £ 196.80

Banking details Bank

Barclays Bank PLC Sort Code 20-78-42 00906662 Account Number

IBAN BIC BUKBGB22
Terms & conditions https://www.viamed.co.uk/terms

GB05BUKB20784200906662

Terms: Net 30 days from date of invoice.

Terms: Net 30 days from date of invoice.
Full invoice amount to be credited to our account net of all bank charges.
Claims: Please claim non delivery within 7 days of invoice.
Shortages or damage within 3 days of receipt.
Claims after these times cannot be entertained.
Title to goods does not pass until payment in full has been received.

Page 1