

COMPANY OPERATING PROCEDURES

Critikon / Sensormedics 12-pin Socket 0010781

VM3/COP/36.32

Date: 3-Jul-02

Revision date: 9-May-11

Issue: 3



Each kit should contain the following components:

Thorkom 12-pin socket		
Qty	Description	Part No.
1	Thorkon 12-pin socket kit	
(1)	 Rear nut	
(1)	 Strain relief	
(1)	 Collett	
(8)	 Sockets	
(1)	 Socket housing	

Procedure:

- 1) Place all the above items into a small plastic bag (9950203).
- 2) Using Lotus Smart Suite – Approach STKLABELS4, print a label as shown below and attach it to the bag and place in the appropriate docket.

0010781
Probe Socket Kit. P890RA.
Critikon/ Sensormedic.