
## anaesthesia

TOF<sup>3D</sup> neuromuscular transmission monitor accessories





5750107

2510103

Overall length: 33.5 cm, each flying lead: 29.5 cm.

Stimulation Cable with pinch clamps\*

Part Number: 2590202

Date: 07/22

## anaesthesia

TOF<sup>3D</sup> neuromuscular transmission monitor accessories





**TOF**<sup>3D</sup> Neuromuscular **Transmission Monitor** Part No.: 2510091





REF Part No. Description 5750110 2520113 Pole Mount

Premium Variant A For use with vertical/horizontal poles and standard medirail. Fits poles up to 4.0 cm diameter or rails up to 5.3 cm.



Acceleration (AMG) Transducer

5750110 Variant B

2520114

Pole Mount Standard

For use with vertical/horizontal poles.



5750116

2530120

Sealing Covers (Pack of 3)

For use with main patient cable, usually used to protect the skin temperature sensor port when not in use. One sealing cover is provided with 2510091.



Cable

When affixing mounts to the unit, push upwards to 'click' securely into position.



Stimulation Skin Temperature

Sensor

5750109

2530121

Sealing Covers (Pack of 3)

For Data Interface Socket. One Sealing Cover included with 2510091.