

Ms Ivana Dordová
ASQA a.s.
Kosor 297
252 26 Kosor
Praha 10, 252 26

Oxygen Hoods

Designed for the short-term stable delivery of oxygen, for use with neonatal and paediatric patients.

Disposable Oxygen Hoods

Consisting of flat sheets of transparent vinyl with built in support frames and a drape neck opening, these hoods are supplied flat packed for compact storage and are easily constructed.


- See through environment provides full patient visibility
- Drape neck opening provides easy access to the patient
- Single patient use means reduced risk of cross-infection
- Flat packed for easy storage and transportation

Standard Hood

The standard oxygen hood has two 22mm inlet ports, a full vinyl bottom and a drape neck opening, and comes in a variety of sizes.

Other style variants of the disposable oxygen hoods are also available, including: No-Bottom variations, NICU, Low Birthweight...

Standard Hood



Please Note:

Oxygen Hoods should only to be used in conjunction with oxygen monitoring.

Enclosed you will find information relating to the Teledyne AX300 Oxygen Analyser and the Teledyne MX300 Oxygen Monitor, available from Viamed.




Please do not hesitate to contact a member of our sales team should you wish to discuss any of these products in further detail.

01535 634542

9:00 – 17:30 Monday to Friday

Yours sincerely



Zoey Teal
Office Administrator

22/5/17

