

VIAMED

Suppliers of Electromedical Equipment

2018



Leeds Liverpool Canal, Kildwick

January

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

February

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	
F	2	9	16	23	
S	3	10	17	24	
S	4	11	18	25	

March

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	31
S	4	11	18	25	

April

M	30	2	9	16	23
T		3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

May

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

June

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	

July

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

August

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	31
S	4	11	18	25	
S	5	12	19	26	

September

M		3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

October

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

November

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	
S	4	11	18	25	

December

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30



VIAMED



15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk

VM-2160 SMARTsat

High Performance

Pulse Oximeter


The very latest and innovative signal processing technologies and algorithms enable precise measurements, even under very difficult physiological conditions



- For spot checking or constant monitoring of SpO₂ and pulse rate
- Now incorporating SMARTsat SpO₂ technology
- Clear and bright, colour TFT display
- Robust, ergonomic design with durable silicone protective cover
- Includes a SMARTsat soft silicone sensor (other sensors also available - see below)
- Stores and restores up to 480 hours of trend data
- Adjustable high and low alarms
- Antibacterial keypad



**SMARTsat
Reusable Adult
Finger Sensor**



**SMARTsat
Reusable Adult
Silicone Soft
Sensor**



**SMARTsat
Reusable
Silicone Wrap
Sensor**



**SMARTsat
Reusable
Ear Sensor**



**SMARTsat
Reusable
Paediatric
Silicone Soft
Sensor**



**SMARTsat
Reusable
Silicone Wrap
Sensor**



VIAMED

Tel: 01535 634 542

Specifications Subject to Change
Part Number: 0196000

