


infant resuscitation



Tom Thumb Infant Resuscitator

This easy to use, compact and robust infant resuscitator, has been designed to enable safe, controlled and fatigue-free resuscitation

- Compact and easy to use
- Complete control over PIP and PEEP
- Safe and consistent resuscitation
- Fatigue-free operation



The Tom Thumb Infant Resuscitator has been designed to replace manual ventilation by the bag and mask method, which can often become tiring and inconsistent. The resuscitation is carried out by using the T-piece occlusion method, thus providing complete control and increased safety over infant resuscitation.

The NeoPEEP neonatal resuscitation circuit, supplied with the Tom Thumb, incorporates a variable PEEP valve and facilitates the delivery of safe and consistent breathing patterns to neonates; by enabling the delivery of consistent PEEP to assist in maintaining functional residual capacity.

Tom Thumb Infant Resuscitator



TT480 Infant Resuscitator



Rear of Tom Thumb with rail clamp

Part Number

0310030

TT480 for connection to an existing flowmeter or air/oxygen blender

Includes:

Variable blow off valve

Fixed blow off valve (predefined, e.g. at 45 cm H₂O)

Integral medirail mounting clamp. Pole mounting clamp also available

NeoPEEP; variable PEEP neonatal resuscitation circuit (single patient use)

Single use silicone round facemask, size 1

Note: This version of the Tom Thumb has a 'Christmas Tree' inlet for the connection to an external flowmeter/blender via oxygen tubing

Accessories

NeoPEEP neonatal resuscitation circuit (single use) - 1m, Box of 20

3210013

NeoPEEP neonatal resuscitation circuit (single use) - 1.5m, Box of 20

3210023

Adapter: 15 mm I.D. - 15mm I.D. - *Reusable*; for connecting the NeoPEEP to the Tom Thumb

Single Unit 0120140

Pack of 10 0120150

Pack of 50 0120152

Adapter: 15 mm I.D. - 15 mm I.D. - *Disposable*; for connecting the NeoPEEP to the Tom Thumb

Single Unit 0120141

Pack of 20 0120151

Box of 100 0120153

Single use silicone round facemask - size 0 - 40 mm (Height) x 27.4 mm (Inner Width) x 48.6 mm (Diameter) - Connector Diameter 15 mm O.D.

Single Unit 3210070

Pack of 20 3210073

Box of 100 3210076

Single use silicone round facemask - size 0.5 - 44.4 mm (Height) x 32 mm (Inner Width) x 54 mm (Diameter) - Connector Diameter 15 mm O.D.

Single Unit 3210067

Pack of 20 3210068

Box of 100 3210069

Single use silicone round facemask - size 1 - 47.9 mm (Height) x 40.5 mm (Inner Width) x 60 mm (Diameter) - Connector Diameter 15 mm O.D.

Single Unit 3210071

Pack of 20 3210074

Box of 100 3210077

Annual service - Tom Thumb

0380000

Service and calibration of instrument to original specification.

Includes all 'O' rings, additional parts are chargeable.

Calibration certificate included.

Other products in our Infant Resuscitation range include: Infant Resuscitation Cabinet, Apgar Timer, Air/Oxygen Blenders

We also supply a range of equipment for use during phototherapy, infant warming and oxygen monitoring.

€ 0086

Specification subject to change

Distributed by:

 **VIAMED**

Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582
Email: info@viamed.co.uk Website: www.viamed.co.uk

