

R-43V Oxygen Sensor Sample Test

See issue #118153

The R-43V Oxygen sensor samples have been tested in air and Oxygen. A calibrated digital volt meter was used to record the outputs (CE206).

Sample 1:

Out of bag test:
26.7mV

2 hours out of bag test:
26.7mV

100% Oxygen test:
131.5mV

Sample 2:

Out of bag test:
27.7mV

2 hours out of bag test:
27.7mV

100% Oxygen test:
131.9mV

Sample 3:

Out of bag test:
26.7mV

2 hours out of bag test:
28.2mV

100% Oxygen test:
134.4mV

Sample 4:

Out of bag test:
32.3mV

2 hours out of bag test:
32.1mV

100% Oxygen test:
152.7mV