

Resuscitation Cabinet Delivery Checklist

Delivery to:

North Manchester General Hospital Delaunay's Road Crumpsall Manchester M8 5RB	Notes: Please notify Dave Gruber (Project Manager - Estates Dept.) on arrival - Tel: 0161 720 4896
---	--

Customer Order No. CAP/004437

Viamed Order No. ORD69490

Part No: 0310002 Description: Infant resuscitation cabinet

Serial No: E348 Item No: 1 of 1

Description	Serial No	Qty	Packed	Rec'd
Cabinet with shelf and rails	E348	1	PC/SH	
Resuscitation bed, 35cm (14") flat	-	1	PC/SH	
Mattress	-	1	PC/SH	
Apgar timer inc. bracket, hex key and batteries	00736	1	PC/SH	
Mounting bracket for APGAR timer	-	1	PC/SH	
Low suction controller, Therapy Equipment 7725	043436	1	PC/SH	
Suction hose, 3m	-	1	PC/SH	
Vacsax suction receiving jar	-	1	PC/SH	
Vacsax disposable liner	-	1	PC/SH	
Vacsax taper	-	1	PC/SH	
Rail clamp with 'V' for mounting Vacsax suction jar	-	1	PC/SH	
Suction tubing, 1m	-	1	PC/SH	
Size 3 plastic storage bin	-	3	PC/SH	
Tom Thumb TT480 inc. patient circuit	0401160	1	PC/SH	
Air/oxygen blender inc. mounting bracket, flowmeter	B1311034	1	PC/SH	
Oxygen hose for blender 3m	-	1	PC/SH	
Air hose for blender 3m	-	1	PC/SH	
Ceratherm 600-3 radiant warmer with mounting arm	468-3-0314	1	PC/SH	
Operator's Manual for cabinet v9.6	-	1	SH	

Please record below any omissions, damage or issues and notify Viamed immediately:

With the exception of any issues noted above, this shipment is accepted complete and in good condition.

Delivered by: (Signature) _____ Print Name: _____

Received by: (Signature) _____ Print Name: _____

Date: _____