

R&D Specifications Sampling Lines Type H for Sedana

Viamed part number	Bluepoint part number	Notes	Material Opaque or clear	Line length	Connectors	Internal diameter	External diameter	External surface, smooth or ridged?	Raised yellow sampling post	For use with anaesthesia	Patent issues	Usage	Packaging
4421515	8090121312V	Same material as 4421514	Opaque	2.5m	LLM - LLM	1.00 mm	2.6 mm	Ridged	No	Requires testing validation	No	Up to 72hrs	Clear plastic pouch with over printed label
4421516	8090121313V	New material	Opaque	2.5m	LLM - LLM	1.00 mm	2.6 mm	Ridged	No	Yes, but still to test and validate	No	Up to 72hrs	Clear plastic pouch with over printed label

To be packed in boxes of 10

For production volumes would be preferential to bulk pack, say 100, 500 or 1000 per box.