



## MaxBlend Lite Compatibility

This leaflet details the characteristics required for oxygen/air blenders to be compatible with the **MaxBlend Lite**.

All medical oxygen air blenders have similar performance requirements in order to comply with international standards (ISO11195 Gas mixers for medical use - Stand-alone gas mixers). However, the physical shape of the blender varies from model to model.

The following diagram detail the blender requirements necessary for compatibility with the **MaxBlend Lite**. If you are unsure if your blender meets these requirements, please contact Viamed for assistance.

