

«TITLE»
«DEPARTMENT»
«COMPANY»
«ADDRESS1»
«ADDRESS2»
«ADDRESS3»
«CITY» «STATE»
«ZIP»

Ref: 1490000/20190806/315P

February 2020

Dear Sir or Madam,

V1000 Foetal Heart Simulator

The V1000 foetal heart simulator is a handheld device, designed for use by technicians and engineers, and for use in any environment where foetal heart monitors need to be tested; before use, when suspected faulty or after routine servicing.

This simulator allows the operation of foetal heart monitors to be verified and the accuracy of displayed foetal heart rates assessed, and can be used to test monitors using standard transducers and gels.

This unique simulator is easy to use and offers objective testing of foetal heart monitors; rather than subjective traditional methods, or cumbersome and expensive test equipment.



The V1000 is to be used with foetal heart monitors and transducers for:

- **Testing prior to use**
- **Fault finding and diagnosis**
- **Soak testing**
- **Simulation and training**
- **Use as part of a PPM procedure**

Please find enclosed the V1000 leaflet and price list overleaf, that we hope you find of interest.

In order to comply with local standards and in response to customer requests, we also offer the following servicing options for the V1000 Foetal Heart Simulator, see overleaf for further information.

- **Certificate of Calibration Only**
- **Service and Certificate of Calibration**
- **Repair Service**

Should you require any further information on this or any of our other products, please do not hesitate to contact a member of our sales team on 01535 634542.

Yours sincerely

Catrin Hollings
Marketing
sales@viamed.co.uk

Please note prices valid at time of printing.

Should you wish to unsubscribe from our mailings
please email unsubscribe@viamed.co.uk

GBP Price List

Part Number	Description	Price (Ex VAT)
1410000	Foetal Heart Simulator V1000	£625.00
0013110	Carrying Case	£33.20
1420000	Wall Mount	£6.85
1420001	Case Protector	£21.00

Servicing Options

- **Certificate of Calibration Only** – Outputs are tested to verify that they are within specification of the settings and are measured in accordance with UKAS traceable calibration standards. The price for this is **£35.00** plus £10.00 return carriage (excluding VAT) – part number 1480001.
- **Service and Certificate of Calibration** – The following are checked and assessed: Outputs settings, LEDs, membrane keys, intensity control, interface cushion and batteries. The outputs are tested to verify that they are within specification at the settings of: 30, 60, 90, 120, 150, 180 and 210 BPM. A functional test is carried out in conjunction with a foetal heart monitor and the unit is left on a soak test. Upon completion a certificate of calibration is issued. The price for this is **£45.00** plus £10.00 return carriage (excluding VAT) – part number 1480000.
- **Repair Service** – This covers the labour on a cost of per hour basis to repair a returned faulty V1000. We will provide a quotation before any work is carried out.

Returns Process

Should you wish to utilise one of these services, please advise us of the serial number/s of the unit/s and we will provide you with a returns reference number. When returning units to us, please may we ask that you ensure that the packing is adequate, batteries are removed and we advise that you use a traceable delivery service.

We recommend, so we can easily recognise your returns, that you use the following formatting for your returns address:

Servicing / Returns Department SRS*****
Viamed Ltd.
15 Station Road
Cross Hills
Keighley
West Yorkshire
BD20 7DT

Please note prices valid at time of printing.

Should you wish to unsubscribe from our mailings
please email unsubscribe@viamed.co.uk