

Department Manager

<<Department>>

«COMPANY»

«ADDRESS1»

«ADDRESS2»

«ADDRESS3»

«CITY»

«ZIP»

Ref: 20200218/536

\*Date\*

Dear Sir or Madam,

### Are you looking for an alternative to the Criterion Oxichheck?

#### The Viamed Oxygen Monitoring Range includes the OxiQuant S and Handi+ Oxygen Analysers

Viamed continually reviews its ranges and has added the **OxiQuant S** and **Handi+ Oxygen Analysers** to its offering.

These units can be used to verify oxygen concentrations in gas mixtures, delivered by: gas supply lines, oxygen concentrators, blenders etc.

The following table highlights the features and benefits of these compact and easy to use units.



#### OxiQuant S

Part Number: 0111275

£155.00

#### Handi+

Part Number: 0111250

£172.00

|                                                                                                                |                        |                              |
|----------------------------------------------------------------------------------------------------------------|------------------------|------------------------------|
| <b>Compact Size</b>                                                                                            | ✓                      | ✓                            |
| <b>Low Cost</b>                                                                                                | ✓                      | ✓                            |
| <b>Sensor</b><br>(replaceable sensor with integrated battery)                                                  | ✓                      | -                            |
| <b>Disposable</b> (integrated sensor and battery)                                                              | -                      | ✓                            |
| <b>Manual Calibration</b> (flexible)                                                                           | ✓                      | -                            |
| <b>Automatic Calibration</b>                                                                                   | -                      | ✓                            |
| <b>Auto Off</b>                                                                                                | 3 minutes              | 80 seconds                   |
| <b>Quick Response Time</b>                                                                                     | ✓                      | ✓                            |
| <b>Warranty</b><br>- <b>Analyser</b><br>- <b>Sensor</b>                                                        | 24 months<br>15 months | 24 months<br>(complete unit) |
| <b>T Adapter</b>                                                                                               | ✓                      | ✓                            |
| <b>Flow Diverter</b>                                                                                           | ✓                      | ✓                            |
| <b>22mm Oxygen Stem Connector</b>                                                                              | ✓                      |                              |
| <b>Oxygen Concentrator Adapter Cap</b>                                                                         |                        | ✓                            |
| <b>Accurate Readings – resolution to 0.1%</b>                                                                  | ✓                      | ✓                            |
| <b>Used to determine oxygen concentrations in gas mixtures, delivered by oxygen concentrators, blenders...</b> | ✓                      | ✓                            |
| <b>Certified for spot checking oxygen percentages during <u>patient treatment</u></b>                          | -                      | ✓                            |
| <b>CE Class</b>                                                                                                | I<br>CE                | II<br>CE0123                 |

## Oxygen Monitoring Catalogue...



The range also includes other oxygen analysers and monitors to measure oxygen levels in a number of different situations and applications, for example in ventilator circuits, oxygen hoods/tents, incubators, anaesthesia machines and CPAP circuits. These products can also be used to test medical equipment and gas supply systems to verify that the correct level of oxygen is being supplied.

Including the **MySign O Oxygen Monitor** which is a highly precise hand-held monitor for clinical, emergency and home ventilation environments. It has a sturdy design which is suitable for paramedics and pre-hospital applications providing immediate access to patients' oxygen concentration parameters, therefore offering safety, even in critical situations.

For further information, please visit our website or for a copy of our latest UK Oxygen Monitoring Catalogue: [www.viamed.co.uk/oxygenmonitoringcatalogue](http://www.viamed.co.uk/oxygenmonitoringcatalogue)

The catalogue also includes the **extended** range of Medical Oxygen Sensors providing the **widest range** and the **latest technology** to offer the best value for money when using OEM and compatible oxygen sensors.

Should you require additional information on applications and compatibility, please visit our website – [www.viamed.co.uk](http://www.viamed.co.uk) or go straight to our cross reference [www.viamed.co.uk/xref](http://www.viamed.co.uk/xref). For any further information please contact a member of our sales team on:

# 01535 634542

(Between the hours of 09:00 and 17:30)

Yours sincerely

Catrin Hollings  
Marketing