

Test Report OOM108 Oxygen Sensor

2

Test Place:	LRI PICU	Tested by:	M. Robinson
Ventilator Master Data			
Ventilator type:	Servo-L	Ventilator Serial Number:	71613
Ventilator Software Version:	V 8.00.01	Ventilator Manufacturing Date:	16022013
Sensor Master Data (currently used / original one)			
Sensor type:	Getinge O2cell	Sensor Part Number:	6640044
Sensor Serial Number:	AA122184	Others:	
Sensor Master Data (tested sensor)			
Sensor type:	OOM108	Sensor Part Number:	E01-00-0075
Sensor Serial Number:	AA600018	Others:	

Test Assignment

Check of oxygen sensor OOM108 in MAQUET ventilator regarding to compatibility, usability, installation, calibration and measurement accuracy.

Checklist

Pos.	Check Item	Yes	No	Comments
1	Sensor usability Sensor usability is same as the original one?	✓		
2	Sensor installation Sensor installation is same as the original one?	✓		
3	Sensor calibration Sensor calibration was successful in the ventilator?	✓		
4	Oxygen Measurement Sensor works without issues in ventilator?	✓		

Final Conclusion:

The oxygen sensor type OOM108 tested in MAQUET ventilator is compatible to the currently used oxygen sensor made by MAQUET.

☒ Yes

☐ No

Date:

15/8/23

Signature:

[Signature]

Test Report OOM108 Oxygen Sensor

①

Test Place:	LRI PICU (w/shop)	Tested by:	MJ Robinson
Ventilator Master Data			
Ventilator type:	Servo-i	Ventilator Serial Number:	18328
Ventilator Software Version:	V8.00.01	Ventilator Manufacturing Date:	(Purchased June 2007)
Sensor Master Data (currently used / original one)			
Sensor type:	O ₂ fuel cell (Geringe)	Sensor Part Number:	06640044
Sensor Serial Number:	AA116453	Others:	
Sensor Master Data (tested sensor)			
Sensor type:	OOM108	Sensor Part Number:	EO1-00-0075
Sensor Serial Number:	AA600018	Others:	

Test Assignment

Check of oxygen sensor OOM108 in MAQUET ventilator regarding to compatibility, usability, installation, calibration and measurement accuracy.

Checklist

Pos.	Check Item	Yes	No	Comments
1	Sensor usability Sensor usability is same as the original one?	✓		
2	Sensor installation Sensor installation is same as the original one?	✓		
3	Sensor calibration Sensor calibration was successful in the ventilator?	✓		
4	Oxygen Measurement Sensor works without issues in ventilator?	✓		

Final Conclusion:

The oxygen sensor type OOM108 tested in MAQUET ventilator is compatible to the currently used oxygen sensor made by MAQUET.

☒ Yes

☐ No

Date:

14/8/23

Signature:



	Revision:	Created - Erstellt:	Page - Seite:
C:\Users\E339512\AppData\Local\Microsoft\Windows\NetCache\Content.Outlook\YEPULFMN\OOM108_compatibility_check-0.docx		MS	1 / 1