



Testimonial Medker EEG Sensor

Medker Tech

User experience report Switzerland

We compared the Medker sensors with the original and quickly found that the values are equivalent.

Thanks to the latest generation chip, the initial signal strength is faster, the handling of the sensor is easy and the adhesive property is strong (3M adhesive).

Prof. Dr. Urs Eichenberger, MD, Head Department of Anaesthesia, Balgrist University Hospital of Zurich, Switzerland

Zurich, 10.05.2023

A handwritten signature in blue ink that appears to read "U. Eichenberger".