

Ms Loe Mair
Manager
Neonatal Unit
Pembury Hospital
Tonbridge Road
Pembury
Tunbridge Wells, Kent TN2 4QJ

March 2012

Dear Ms Mair,

POSEY Pulse Oximetry Sensor Wraps



The Posey pulse oximetry sensor wraps replace adhesive tape to secure a wide variety of flexible pulse oximetry sensors. The wraps hold the sensor firmly to prevent misalignment of optical components, whilst blocking out excessive ambient light.

The soft foam tape-less fastener does not stick to the skin, allowing for easy application and removal; this also helps to preserve delicate skin. The model ref. 6554 is a universal design and may be cut to length, the model ref. 6554P (premie) has been specifically designed for neonatal applications.

Part Number	Description	Price	Pack Qty
0021013	Model: 6554 Application: Neonatal Foot, Infant Toe or Adult/Paediatric Finger Pack of 12	£6.00	1
0021014	Model: 6554 Application: Neonatal Foot, Infant Toe or Adult/Paediatric Finger One Case (48 boxes of 0021013)	£288.00	48
0021018	Model: 6554P Application: Neonatal Foot – Premie Pack of 12	£8.00	1
0021019	Model: 6554P Application: Neonatal Foot - Premie One Case (48 boxes of 0021018)	£384.00	48

I hope you find the above of interest; should you require a free of charge sample or further information, please do not hesitate to contact a member of our sales team on 01535 634 542 between 09:00 and 17:30.

Yours sincerely



Steve Hardaker
UK Sales Manager