

Catrin Hird

From: Catrin Hird <catrin.hird@viamed.co.uk>
Sent: 15 March 2012 13:45
To: Catrin Hird
Subject: Viamed - VM-2500 Capnograph with Multigas

Dear Catrin Hird,

High performance and versatile anaesthetic agent monitoring

The ultra-compact, light and easy to handle design of the VM-2500-MG Multigas Monitor is the perfect flexible and mobile monitor for identifying and quantifying the five most important anaesthetic agents as well as other gases and parameters.

- **Innovative micro-optic technology**
- **Direct mainstream measurement without time delay**
- Reaches **full-scale accuracy** within **1 minute** to full specification
- **No service or calibration required**; reduces running costs and down time.
- **Bright and clear organic LED display**

Utilising innovative technologies and state-of-the-art design, this CO₂ and SpO₂ monitor provides extremely accurate and reliable measurement for intubated or non-intubated patients.

The VM-2500-MG capnograph utilises the advanced IRMA™ AX+ Analyser; to precisely determine gas concentration in the mixture – with state-of-the-art NDIR technology with up to 9-channel gas type analysis in the 4 – 10 µm range. The VM-2500-MG is designed to overcome the shortfalls of conventional sidestream technologies; namely water and secretion handling, calibration and service costs. The IRMA™ AX+ Analyser is factory calibrated and requires no integration procedures or associated expenses, just 'PLUG-IN and MEASURE...™'.

Measuring a range of vital parameters; FiCO₂, EtCO₂ Respiration Rate, SpO₂ and Pulse Rate as well as the five most important anaesthetic agents - Halothane, Isoflurane, Enflurane, Sevoflurane and Desflurane. With no maintenance or calibration required, this high performance capnograph is a cost effective system for patient monitoring during anaesthesia, even in the most difficult situations.

I have attached a copy of the Capnograph with Multigas leaflet for your reference, should you require any further information please do not hesitate in contacting me.

Yours sincerely

Kind regards

Catrin Hird
Marketing Co-ordinator
Viamed Ltd.

<http://www.viamed.co.uk>
Email catrin.hird@viamed.co.uk
Tel: +44 (0)1535 634 542
Fax: +44 (0)1535 635 582

Skype: viamed.catrin.hird

Twitter: twitter.com/ViamedLtd
Facebook Page: Search for Viamed Ltd

Company registered in England, No. 01291765.

E-mail Disclaimer
