

Catrin Hird

From: Catrin Hird <catrin.hird@viamed.co.uk>
Sent: 19 March 2012 11:40
To: Catrin Hird
Subject: Viamed - NeoPEEP

Dear Catrin Hird,

NeoPEEP - Variable PEEP Neonatal Resuscitation Circuit

The NeoPEEP neonatal resuscitation circuit incorporates a variable PEEP valve and facilitates the delivery of safe and consistent breath patterns to neonates; by enabling delivery of consistent PEEP to assist in maintaining functional residual capacity.

This resuscitation circuit will connect directly to any resuscitator with a standard 15m I.D. patient outlet, such as the Hill-Rom/Draeger resuscitation units, as well as standard endotracheal tubing. We also have connectors available to adapt the circuit to fit most other resuscitation units which feature a non-standard patient outlet, such as our Tom Thumb Infant Resuscitator and the Fisher & Paykel NeoPuff.

For your reference I have attached a leaflet which includes further information. Should you require any additional information please do not hesitate in contacting me.

Yours sincerely

Kind regards

Catrin Hird
Marketing Co-ordinator
Viamed Ltd.

<http://www.viamed.co.uk>
Email catrin.hird@viamed.co.uk
Tel: +44 (0)1535 634 542
Fax: +44 (0)1535 635 582

Skype: viamed.catrin.hird

Twitter: twitter.com/ViamedLtd
Facebook Page: Search for Viamed Ltd

Company registered in England, No. 01291765.

E-mail Disclaimer

The information in this email is confidential and may be legally privileged.
It is intended solely for the addressee. Access to this email by anyone else is unauthorized.
If you are not the intended recipient, any disclosure, copying, distribution or any action taken or omitted to be taken in reliance on it, is prohibited and may be unlawful.

If you feel you have received this email in error or would like to be removed from our mailing list, please reply with delete in the subject line or call +44 (0)1535 634542
