

VIAMED

Suppliers of Electromedical Equipment

2013



East Riddlesden Hall, Riddlesden

January

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

February

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	
S	2	9	16	23	
S	3	10	17	24	

March

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31

April

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

May

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	31
S	4	11	18	25	
S	5	12	19	26	

June

M		3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

July

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

August

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	31
S	4	11	18	25	

September

M	30	2	9	16	23
T		3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

October

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

November

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31

December

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29



VIAMED



15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk

Capnography

Viamed offer the VM-2500 range of high performance Capnographs addresses the growing emphasis on the use of Capnography during resuscitation and intubation, inline with recommendations of NAP4 report



VM-2500-M
Part Number: 4410500
EAN-13: 5051826002711

VM-2500-MG Multigas/SpO₂ Monitor

High Performance Capnograph with either Mainstream or Sidestream Measurement

- Fast warm-up time of <10 seconds to full scale accuracy
- Portable Capnograph including Pulse Oximetry as standard.
- Mains and battery powered, rechargeable Li-Poly or standard AA alkaline batteries provided.
- No service or calibration required; reduces running costs and down time.
- Bright and clear organic LED display; can be read from any angle.
- 400 hours trending data, downloadable to PC. Includes PC software as standard.

VM-2500-S
Part Number: 4410520
EAN-13: 5051826002728



VM-2500-MG
Part Number: 4410550
EAN-13: 5051826003190

VM-2500 CO₂/SpO₂ Monitor

High performance and versatile anaesthetic agent monitoring

- Innovative micro-optic technology
- Direct mainstream measurement without time delay
- Warm-up time: < 1 minute to full specification
- Maintenance and calibration-free technology
- Colour OLED display
- Application range from neonates to adults



Short Term (Single Patient Use)

Nasal Sampling Line For triage and emergency services applications
Adult 4420560 Box of 40 Length: 2m

Long Term (Single Patient Use)

Nasal Sampling Cannula Length: 2m

Adult 4420710 Box of 25

Paediatric 4420711 Box of 25

Infant 4420713 Box of 25

Nasal Sampling Cannula with O₂ Delivery Length: 2m

Adult 4420720 Box of 25

Paediatric 4420721 Box of 25

Nasal/Oral Sampling Cannula Length: 2m


Adult 4420730 Box of 25

Paediatric 4420731 Box of 25

Nasal/Oral Sampling Cannula with O₂ Delivery Length: 2m

Adult 4420740 Box of 25

Paediatric 4420741 Box of 25



Alternative sensors available:




Adult Soft Silicone Sensor*




Paediatric Soft Silicone Sensor



Finger Sensor



Wrap Sensor



Ear Sensor

Range of Disposable Sensors also available

* Standard Sensor Provided



Colour OLED Display



ISO 9919:2009



Clinical Validation



VIAMED

Tel: 01535 634 542

Specifications Subject to Change
Part Number: 0196000

