



## SpiroTrue Flowsensors



For flow measurement of gas delivered by anaesthesia systems and ventilators

### SpiroTrue A

**Hot-wire flowsensor for use with Dräger ventilators**

- Compatible to the Dräger Spirolog flow sensor: Part no. 8403735.
- Hot-wire flow sensor for use with Dräger equipment: Dräger Spirolog, Anemone, PM8030, PM8060, Eva, Cierco, Cato, Julian, Evita 4, Evita XL, Primus...
- Supplied in a pack of five, packed individually in sterilizable packaging.
- The sensor can be sterilised with ETO before usage.
- Several disinfection cycles with 70% Ethanol solution are possible. The sensor can be used as long as equipment calibration is possible.

Viamed Part Number: 4310001



## Autoclavable

### SpiroTrue A (PC)

- Compatible to the Dräger SpiroLife flow sensor: Part no. MK 01900.
- Hot-wire flow sensor for use with Dräger equipment: Dräger Spirolog, Anemone, PM8030, PM8060, Eva, Cierco, Cato, Julian, Evita 4, Evita XL, Primus...
- Supplied in a pack of five, packed individually in sterilizable packaging.
- The sensor can be sterilised with ETO before usage.
- Several disinfection cycles with 70% Ethanol solution are possible. The sensor can be used as long as equipment calibration is possible.


Viamed Part Number: 4310002

### SpiroTrue H

**Differential pressure flowsensor for use with Hamilton ventilators**

- Replacement for Hamilton H7400 Accu Tach flow sensor.
- For use with Hamilton H7200 VR1, Galilio, Amadeus, Veolar, Raphael...
- SpiroTrue H is supplied in a pack of six, packaged individually in sterilizable packaging.

Viamed Part Number: 4310003



# VIAMED

**Viamed Limited** · 15 Station Road · Cross Hills  
Keighley · West Yorkshire · BD20 7DT · United Kingdom  
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)



???????

Specifications subject to change

Date: 07/13



BS EN ISO 9001:2008  
ISO 13485:2003