



NeoMask Phototherapy Eye Masks

The soft, safe and secure solution for eye protection during neonatal phototherapy

- Adjustable band for optimum comfort and a secure fit
- Easy to use design
- Ergonomically designed eye protection pad to reduce pressure and increase light protection
- Ultrasonically welded tabs - no adhesives for a secure fastening



Part Number: 1114015

Regular


?? - ?? cm occipital - frontal circumference



Part Number: 1114016

Small

?? - ?? cm occipital - frontal circumference



Part Number: 1114017

Extra Small

?? - ?? cm occipital - frontal circumference



The EyeMax2 has been designed to assure optimum comfort and eye protection for neonates during phototherapy treatment. The single-piece, wrap-around design is secured at two points on either side of the neonate's head and can be independently adjusted to prevent unwanted movement and allows a comfortable and secure fit.

The large eye-protection pad has been ergonomically designed to reduce pressure and increase light protection. The headband; however, allows blue light to penetrate, so that the head continues to receive phototherapy treatment.

All phototherapy eye masks from Viamed are supplied in packs of 20

Other phototherapy products available from Viamed

An alternative to the traditional phototherapy eye mask is the phototherapy light shield. These light shields are manufactured by Viamed using Amber 300™ material, and are available in a variety of standard designs and sizes. Light shields can also be designed to meet exact specifications, should the standard range not meet individual requirements.

Distributed by:



VIAMED

Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582
Email: info@viamed.co.uk Website: www.viamed.co.uk



BS EN ISO 9001:2008
ISO 13485:2003



Part Number: 1114015
Date: 01/01/00

FS 28344 MD 78787