


pulse oximetry



MD300-C

Finger Pulse Oximeter

Ultra-compact, lightweight pulse oximeter with LED display

- Low cost
- Easy to use; 1-button operation
- Automatic power-off when finger is removed

**Economical
pulse oximeter
for general
applications**



Actual Size

The MD300-C can be used in a variety of applications, from use in hospitals, paramedic ambulances, clinics and dental surgeries, through to spot-checking before, during and after exercise.



VIAMED

MD300-C Finger Pulse Oximeter

Measuring Range

| | |
|------------------|--------------|
| SpO ₂ | 70 - 100% |
| Pulse Rate | 30 - 250 bpm |

Accuracy

| | |
|--|------------|
| SpO ₂ (70 - 100%) | ± 2 digits |
| SpO ₂ Resolution | 1% |
| Pulse Rate; beats per minute (30 - 99 bpm) | ± 2 bpm |
| Pulse Rate; beats per minute (100 - 250 bpm) | ± 2% |
| Pulse Rate Resolution | 1 bpm |

Display

| | |
|----------|-------------|
| Features | LED display |
|----------|-------------|

| | |
|----------------|---|
| Data displayed | Oxygen saturation, pulse rate, pulse indicator, bar graph |
|----------------|---|

| | |
|------------|---------------------------------|
| Indicators | Signal strength, battery status |
|------------|---------------------------------|

Environmental Conditions

| | |
|-----------------------|----------------|
| Operating temperature | 5°C to 40°C |
| Storage temperature | -20°C to +55°C |

Regulatory Approvals

| | |
|-----------|--------------|
| Standards | CE Class IIa |
|-----------|--------------|

Assembly

| | |
|--------------|---|
| Dimensions | 56 mm (L) x 32 mm (W) x 30 mm (H) |
| Weight | Approx. 50 g (including batteries) |
| Power Supply | 2 x AAA/LR03/MN2400 1.5V alkaline batteries |
| Battery Life | Approximately 16 hours of continuous operation. Automatically powers off when finger is removed for longer than 8 seconds. |

Part Numbers

| | |
|---|---------|
| MD300-C | 2810010 |
| EAN-13 (GTIN-13) Barcode: 5051826003046 | |
| Standard accessories included: | |
| Lanyard | |
| 2 x AAA Batteries | |
| Operator Manual | |

Optional Accessories

| | |
|--|---------|
| Soft Carry Case, with belt loop, colour - black | 0022190 |
| Soft Carry Case, with belt loop, colour - orange | 0022191 |
| Soft Carry Case, with belt loop, colour - green | 0022192 |
| Soft Carry Case, with belt loop, colour - blue | 0022193 |

Clinical trials carried out in accordance with BS EN ISO 9919:2005 standards.
Materials tested in accordance with ISO 10993-, -5, -10.



Optional carrying cases with belt loop and D-ring

CE 0123

Specifications Subject to Change

Distributed by:



VIAMED

Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582
Email: info@viamed.co.uk Website: www.viamed.co.uk



Part Number: 2890099
Date: 09/14



ISO 13485:2003