

Existing Capnograph customers

ACCEPTED.....
DATE 10.3.15
INITIAL 8

Supplies Manager
Supplies Department
East Lancashire Hospitals NHS Trust
Royal Blackburn Hospital (Trust HQ)
Haslingden Road

Blackburn Lancashire
BB2 3HH

ACCEPTED.....
DATE 6/3/15
INITIAL SH

*will not
be original
docs*

February 2015

Dear Sir or Madam,

**Introducing VersaStream CO₂ Sampling Lines
for use with Viamed VM-2500 Series Capnographs**

Viamed are excited to introduce VersaStream CO₂ sampling lines as an innovative solution for sidestream CO₂ monitoring, offering a one-piece sampling solution with increased performance and low cost. These can be used as an alternative option instead of the original Nomo adapter and Nomolines as the whole line is single-patient use and has an integral, compact filter.

VersaStream CO₂ sampling lines use the latest material technologies to allow the complete removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

Why VersaStream CO₂ sampling lines?

- Cost-effective alternative to Nafion[®] tubes
- The gas-selective non-porous material facilitates fast and complete removal of condensation from the sampling line
- The integrated bacterial filter prevents contamination
- Reliable measured gas values, due to negligible dead space volume of the filter - even in cases of high respiration rates
- Integrated bacteria filter and water removal reduces the number of patient circuit parts and speeds up patient monitoring preparation time

Available in a number of variants including: single nasal sampling line, nasal cannula, nasal/oral cannula and a one piece airway adapter sampling line.

Please find enclosed a leaflet outlining the range along with a price list. Should you have any further questions regarding this, or would like to obtain samples for evaluation, please do not hesitate in contacting our sales team on (01535) 634542 between 9:00am and 5:30pm Monday - Friday.

Yours sincerely

C Hollings
Catrin Hollings
Sales & Marketing Co-ordinator

ATI customers

April 2015 → Ambulance Trusts

PICS

ACCEPTED.....
DATE 10/3/15.....
INITIAL *su*.....

PICS

Supplies Manager
Supplies Department
East Lancashire Hospitals NHS Trust
Royal Blackburn Hospital (Trust HQ)
Haslingden Road

Blackburn Lancashire
BB2 3HH

ACCEPTED.....
DATE 10-3-15.....
INITIAL *su*.....

March 2015

Dear Sir or Madam,

**Introducing VersaStream Gas Sampling Lines for use with
Capnographs and Multigas Monitors with Luer Lock Male or Respironics/Philips connector**

Viamed are excited to introduce VersaStream gas sampling lines as an innovative solution for sidestream gas monitoring, offering a one-piece sampling line solution with increased performance and low cost.

VersaStream gas sampling lines use the latest material technologies to allow the complete removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

Why VersaStream gas sampling lines?

- Cost-effective alternative to Nafion® tubes
- The gas-selective non-porous material facilitates fast and complete removal of condensation from the sampling line
- The integrated bacterial filter prevents contamination
- Can reduce the costs of regular water trap exchange (if present)
- Reliable measured gas values, due to negligible dead space volume of the filter - even in cases of high respiration rates
- Integrated bacteria filter and water removal reduces the number of patient circuit parts and speeds up patient monitoring preparation time
- VersaStream Types available:
 - Luer Lock Male gas sampling line for use with many common types of gas exchange monitors, either low-flow or high flow
 - Respironics/Philips CO₂ sampling line for use with LoFlo™ sidestream modules

Both ranges are available in a number of variants including: single nasal sampling line, nasal cannula, nasal/oral cannula and a one piece airway adapter sampling line.

Please find enclosed a leaflet outlining the range along with a price list. Should you have any further questions regarding this, or would like to obtain samples for evaluation, please do not hesitate in contacting our sales team on (01535) 634542 between 9:00am and 5:30pm Monday - Friday.

Yours sincerely

C Hollings

Catrin Hollings
Sales & Marketing Co-ordinator

Ambulance Trusts : new potential customers

Supplies Manager
Supplies Department
East Lancashire Hospitals NHS Trust
Royal Blackburn Hospital (Trust HQ)
Haslingden Road

Blackburn Lancashire
BB2 3HH

ACCEPTED.....
DATE 10.7.15
INITIAL BW

ACCEPTED.....
DATE 6/3/15
INITIAL SH

February 2015

Dear Sir or Madam,

The increased emphasis in Capnography within Resuscitation Situations

As you may be aware, the Resuscitation Guidelines released by the Resuscitation Council (UK) in October 2010 recommends the use of Capnography in addition to clinical assessment as the most reliable method to confirm placement of a tracheal tube.

To help with this Viamed offer the VM-2500 range, which includes both mainstream and sidestream units. These units are of high specification without the high price tag seen by other units on the market.

Key features of the VM-2500 Capnography range from Viamed include:

- Fast warm-up time less than 10 seconds to full specification
- Rapid response time
- Extremely accurate and reliable measurement
- Clear and bright, colour display offers superior clarity
- For intubated and non-intubated patients
- Range of different SpO₂ sensors available including adult and paediatric soft sensor, standard finger sensor, silicone wrap sensor, range of disposable sensors and ear sensors.



Introducing VersaStream CO₂ Sampling Lines for use with Viamed VM-2500 Series Capnographs

Viamed are excited to introduce VersaStream CO₂ sampling lines as an innovative solution for sidestream CO₂ monitoring, offering a one-piece sampling solution with increased performance and low cost.

VersaStream CO₂ sampling lines use the latest material technologies to allow the complete removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

Why VersaStream CO₂ sampling lines?

- Cost-effective alternative to Nafion® tubes
- The gas-selective non-porous material facilitates fast and complete removal of condensation from the sampling line
- The integrated bacterial filter prevents contamination