

oxygen monitoring mounting options

for use with **TELEDYNE**
ANALYTICAL
INSTRUMENTS

AX300 and MX300 oxygen monitors




AX300

- 0111230** Standard Unit
3.5 digit display, analogue output
- 0111232** 3.5 digit display, RS232 output
- 0111233** 2.5 digit display, analogue output
- 0111234** 2.5 digit display, RS232 output




MX300

- 0111235** Standard Unit
3.5 digit display, analogue output
- 0111237** 3.5 digit display, RS232 output
- 0111238** 2.5 digit display, analogue output
- 0111239** 2.5 digit display, RS232 output



0121213 Mounting Plate for use with the AX300 and MX300



0121180 V Wall Mount




0121181 V Mount Pole Clamp - Small
Clamping Width: 14 - 25mm



0121182 V Mount Pole Clamp - Large
Clamping Width: 16 - 40mm



0121184 V Mount Rail Clamp
Attaches to most types of 30mm (H) x 10mm (D) Medirail



0121200 Universal Mounting Clamp
Attaches to round or flat objects up to 3.8 cm thick. Facilitates positioning of the device at any angle



Connects direct to the monitor



Part Number: ???????
Date: 10/15