

ACCEPTED.....
DATE 4/5/16.....
INITIAL *AK*.....

ACCEPTED.....
DATE 4/5/16.....
INITIAL *AK*.....

Visit us at
Stand K20 at
The Emergency
Services Show

VM-2500 range of high performance Capnographs
address the growing emphasis
on the use of Capnography
during resuscitation and intubation

Start-up
full scale
accuracy
within 10
seconds

VersaStream®
Gas Sampling Lines
Versatile, lower cost, innovative
solutions for sidestream monitoring

Available as Mainstream
and Sidestream units

VIAMED
Tel: 01535 634542
E-mail: EST@viamed.co.uk
www.viamed.co.uk

