


Specification Datasheet

Viamed

SpiroTrue H

OEM Disposable Single Patient Adult Flow Sensor



CE

Part numbers	4310003 pack of 6 4310013 OEM pack of 20
Model	SpiroTrue H – OEM disposable single patient adult flow sensor + tubing set
Primary packaging	PE pouch (bag) with 20 sets includes flow sensor with double PVC (DEHP free) tube connected plus adapter and hook

Technical Data

Measuring principle	Differential pressure measurement
Reproducibility	≤ 10% at flow range 0 to 100l/min
Designed in accordance with	93/42/EEC (Annex I and Annex II)
Housing material	ABS (Acrylnitril-Butadien-Styrol)
Weight	≤ 9.5g (without tubing)
Dead space	≤ 10ml
Connectors	15mm O.D./15mm O.D. in accordance with ISO 5356-1
Conforms to	EN ISO 10993-5
Packaging	Non-sterile
Recommended storage temperature	-20°C to 50°C
Humidity (R.H.)	5% to 95%
Storage period	≤ 3 years (unopened, away from direct sunlight)
Operating temperature	10°C to 40°C
Relative humidity	5 to 95%, non-condensing

Notes: Before using the sensor consult the appropriate ventilator user instructions for more detailed information.

SpiroTrue H is compatible with Hamilton single use flow sensor part number 279331